

HENRY SCHEIN®

OneStep™ + Pro Hb Analyzer

The Henry Schein OneStep+ Pro Hb Analyzer provides an accurate hemoglobin and hematocrit measurement.

1
test

2
results

3
seconds

- Reagent-free microcuvettes
- Microcuvette shelf life of up to 2 years, even after the canister has been opened
- Cleared for use on patients 6 months or older
- Measurement range 5.0 to 25.6 g/dL
- Small 15 µl sample volume
- 3.25" Color LCD display
- CLIA Waived



Reimbursement

Test	Hemoglobin
CPT Code	85018-QW
Blood count; hemoglobin (Hgb)	



HENRY SCHEIN®



Simple 3 Step Procedure



• Performance



• Parameters

Specifications	OneStep+ Pro Hb
Sample Type	Capillary Blood/Venous Blood
Reaction Time	3 Seconds
Measurement Range	5.0 to 25.6 g/dL
Sample Volume	15 µl
Microcuvettes Stability	24 Months
Storage Condition	15 - 35°C / 59 - 95°F
Tests/box	200

• Analyzer

Specifications	OneStep+ Pro Hb
Dimensions	4.19 x 5.94 x 1.51 in
Weight	0.5 lbs.
Interface	3.25-inch Screen External Thermal Printer
Memory	1,000 Results and QC Results
Power	4 AA Batteries*

• Order Information

Products	Item Number
Pro Hb Analyzer	570-1311
Pro Hb Microcuvettes	570-1312
Pro Hb Controls (Level 1,2,3)	570-1313
Pro Hb Optical System Check	570-1322
Pro Hb Starter Kit**	570-1314

* Batteries and USB Power Cable Included.

** Starter Kit Includes 1 Analyzer, 1 Optical System Check, 100 Microcuvettes.

PS-1545 Rev02 2019/02



Contact your Henry Schein Sales Consultant for product information.

1-800-772-4346

www.henryschein.com/medical