

MANUFACTURED BY:



Packaged in metal storage case!



- Volume Control
- Battery Door
- Microphone

Introducing the PHD, the Physician's Hearing Device!


A professional assistive listening device manufactured to FDA medical device quality standards. Your hearing impaired patients will appreciate the value added service of in-office amplification to enhance their hearing during routine office visits.

The PHD is preprogrammed to provide a smooth, broad-band, sound quality to offer assistance to patients who have a mild to moderate high frequency hearing loss.

The PHD serves three major purposes:

- Immediate amplification to facilitate communication between the physician and patient.
- Fulfills current directives with the ADA mandate for effective communication.
- Demonstrates benefits of binaural amplification.

The PHD is specifically designed for clinical use, simple to use and maintain. Offering durable, medical quality components housed in a polypropylene case. This material is hypoallergenic and does not support bacteriological growth. The chrome plated headset and soft ear tips quickly clean with alcohol wipes recommended in clinical practices.



Powered by a standard #13 size hearing aid battery

PHD SPECIFICATIONS

- Battery Life525 Hours
- Battery Drain.....0.59 mA
- Output.....110 dB
- Gain30 dB HF avg.



©2016 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. Note: Participating in a promotional discount program (e.g., points, discount redemptions or other special awards) is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the program requirements. For HS promotions: This promotion offers a bundled discount. The cost of any "free" (or "no charge") items is a discount applied to your other purchases on the same invoice. The reportable value (unit price) of the "no-charge" item will be identified on the invoice. To calculate the applicable discount on the other items on the same invoice, divide the total value (unit price(s) of any no-charge items(s) by the total value (extension price(s) of all other items purchased on the same invoice, and apply the resulting percentage equally to each item to determine net prices. For manufacturer promotions: The prices stated in the invoice or statement may reflect a discount or be subject to a rebate. You must fully and accurately report this stated discount price, or if applicable, any net pricing, after giving effect to any rebates, to Medicare, Medicaid, Tricare and any other federal or State program upon request by any such program. Accordingly, you should retain these records. It is your responsibility to review any agreements or other documents applicable to these prices to determine if they are subject to a rebate.

Available through



To Order: 1-800-372-4346 8am-9pm (et)
To Fax: 1-800-732-7023 24 Hours

ITEM CODE #1286112