

Professional Vacuum Solutions

Dominant 50 — Powerful and reliable

Maximum suction power with Dominant 50

- Fast vacuum build-up and high suction capacity
- Precise vacuum
- Vibration free, quiet functioning
- Appropriate for large volume fluid collection

Henry Schein Mobile#: 240-0019

Henry Schein Portable#: 240-0020



Basic 30 — Versatile and highly efficient

The peak of technology

- Rapid vacuum build-up
- Powerful suction and precise vacuum

Great flexibility

- As a portable model or on a mobile stand with foot on/off switch
- Ideal size – fits into endoscopy tower

Henry Schein Mobile#: 240-0023

Henry Schein Portable#: 240-0024



Vario — Mobile and flexible handling

State-of-the-art technology

- Vacuum production with innovative QuatroFlex™ technology
- Precise vacuum
- Small, light and quiet

Model variety

- Different vacuum versions (Vario 18 AC, Vario 18 AC/DC, Vario 18 AC/DC c/i)
- Models available with integrated battery which can be charged during operation
- 3 different warning signals during battery use

Henry Schein AC#: 240-0021

Henry Schein AC/DC#: 240-0022



Fluid Management System

Enhanced safety and simplicity

- 1500cc and 2500cc liner with handles
- Color coded canisters and liners
- Canister with locking mechanism prevents spills when removing liners
- Inherently allows for safer handling of potentially infectious fluids
- Reusable autoclavable canister options are available



Pump Quick Facts Tool

✘ = ideal
● = possible

| | Medical Suction Systems | | | | | | Application Area | | | | |
|---|--------------------------|--------|-------------------|-----------------|-------------|--------|------------------|---------|---------------|---------------------------|--------------------------|
| | Suction capacity (l/min) | (mmHg) | Max. vacuum (kPa) | Body contouring | Liposuction | Breast | Face, head, neck | Surgery | Room/ward/ASC | Practice/physician office | OR/sterile procedure rm. |
| Professional Vacuum Mobile Suction | | | | | | | | | | | |
| Vario 18 AC | 18 | -563 | -75 | | | | ✘ | | ✘ | ✘ | ✘ |
| Vario 18 AC/DC | 18 | -563 | -75 | | | | ✘ | | ✘ | ✘ | ✘ |
| Basic 30 | 30 | -638 | -85 | | | ✘ | | ✘ | ● | ✘ | ✘ |
| Dominant 50 | 50 | -675 | -90 | ✘ | ✘ | ● | | ✘ | | ● | ✘ |

Call your Henry Schein representative for more information.

Visit www.henryschein.com/medela

For a complementary in office demonstration, visit:
www.henryschein.com/requestinfo