

## German Biopsy Punch Forceps

German Precision Instruments



- Exact bite for clean tissue extraction
- Firm hold of biopsy sample for confident procedures
- Spring grip for optimal handling
- Swivel lock to ensure sample retention
- Popular 10" length
- Outstanding value, one-year unconditional warranty

Made in Germany from high-grade surgical stainless steel, Premier's Biopsy Punch Forceps offer optimal performance and ultimate quality. Through superior craftsmanship, these four classic shapes afford the highest level of clinical care.

See reverse for full product selection.


## German Biopsy Punch Forceps


German Precision Instruments



**9086420 Tischler**


This clinical standard features a unique bite mechanism which grasps tissue and holds it securely during excision.


Bite: 7mm x 3mm 



**9086410 Burke**


This small-toothed forceps allows for sample retrieval with minimal tissue penetration.


Bite: 5mm x 3mm 



**9086430 Baby Tischler**


A mini or Baby Tischler provides a smaller sample. Its rounded bite minimizes patient discomfort and reduces the occurrence of bleeding.

Bite: 4.2mm x 2.3mm 



**9086440 Kevorkian**

This classic pattern offers the largest available jaw and the serrated teeth afford strong grasping ability.

Bite: 9mm x 3mm 

## Monseľ's

Ferric Substrate Solution

Pre-Thickened, Ready-to-use!

**9045055** 12 single application 8ml bottles and 12 rayon tip applicators



[premusa.com/promise](http://premusa.com/promise)