

ENT

Glass Bead Mirror Warmer

Warms Mirrors Prior to Laryngeal Examinations to Prevent Fogging



Affording the ENT Clinician improved and efficient laryngeal examinations, while offering the following benefits:

- **Rapid warm-up**¹ – Ready for use with the day's first patient
- **Thermostatically controlled temperature** – Turn on in the morning and do nothing until leaving for the evening
- **Insulated thermal walls** – Housing of the unit remains cool and safe to touch
- **Low profile, space-saving design** – Less clutter in the examination room
- **Large bead well** – Less spillage, improved stability
- **Hinged lid** – Convenient access to bead well, limits exposure to beads when not in use
- **UL® listed** – Meets or surpasses most institutional safety standards
- **Designed to fit in your existing cabinet.**² – Contact your examination cabinet supplier or surgical supply distributor for details.

Simply immerse a sterile mirror in the hot beads for three to five seconds prior to examination.

Cautions:

- Allow mirror sufficient time to cool before use
- Soiled or contaminated beads must be changed

- 1006131 Mirror Warmer (110V) and Glass Beads
- 1006132 Mirror Warmer³ (220V) and Glass Beads
- 1006105 Replacement Beads

premier[®]

Premier® Medical Products • 888-670-6100 • www.premusa.com