



## Uni-Punch® Seamless Disposable Biopsy Punch

- **Seamless Blade**  
Sharpened on both inside and outside
- **Ergonomic Handle Design**  
Lightweight, stiff plastic material, ribbed handle for improved control
- **Sterile Product**  
Ethylene oxide sterilization (ETO)
- **Mylar and Tyvek Sterile Peel-Pouch**  
Individually wrapped
- **Capped Tip for Added Safety**  
During handling and disposal
- **Packaging Convenience**  
Boxes of 25 and 50, Assorted package



# Uni-Punch®

## Seamless Disposable Biopsy Punch



Nine Different Sizes!  
The largest number available!



(Assorted - 5 each: 2mm, 3mm, 3.5mm, 4mm, 5mm)

### Box of 25

Size (in mm)	Item #
1.5	9033509
2	9033502
2.5	9033511
3	9033503
3.5	9033510
4	9033504
5	9033505
6	9033506
8	9033508
Asst.	9033525

### 1.5mm



### 2mm



### 2.5mm



### 3mm



### 3.5mm



### 4mm



### 5mm



### 6mm



### 8mm



### Box of 50

Size (in mm)	Item #
1.5	9033519
2	9033512
2.5	9033521
3	9033513
3.5	9033520
4	9033514
5	9033515
6	9033516
8	9033518

Available from:

