



Pro Plus 750HS Series Clockface Aneroids

A large face aneroid available in 2 configurations that offer the superb precision, styling, and dependability.

The Pro Plus 750HS series features:

- Large, easy-to-read 6" luminescent dial for high visibility in virtually any light
- 300mmHg CE no-pin stop movement withstands the heaviest use
- Inflation system includes black nylon cuff, latex free inflation bladder, deluxe bulb and valve, and eight foot length (when extended) of coiled latex free tubing
- Contemporary Euro-style ivory casing compliments any decor
- Complete with all necessary hardware and assembly instructions
- Available in Wall or Mobile configurations
- Assembled, inspected, and packaged in the U.S.A. from foreign components

The Pro Plus 750W Wall Aneroid includes:

- Wall bracket that swivels through 110° rotation for easy viewing with integral cuff storage

The Pro Plus 752M Mobile Aneroid includes:

- Height adjustable Spider-Leg™ caster stand with range of 45" to 54" (from bottom of base to top of unit)
- Integral ABS cuff storage compartment
- Manometer angled back 15° for easy viewing
- 18" base diameter

**Some assembly required*

PRIVATE LABELING	
COMPONENT	PACKAGING
Gauge & Cuff	Label



Feature	Benefit
Nissei Made Aneroid Gauge EXCLUSIVELY - The MOST experienced and respected Asian Gauge manufacturer	Unsurpassed accuracy, quality, and reliability. Our gauges are randomly tested to 30,000 cycles, 3x MORE than required under the National Standard (ANSI SP 10)
100% visually inspected for 0 point, 15% re-inspected for accuracy.	To ensure conformance with National Standard ANSI SP 10
Gauge features dampening mechanism	Reduces "needle bounce" and promotes longer life.
Gauge features Luminescent dial	Easier to read in low light conditions
Gauge crystal specially designed to seal out dust	Prolongs life of movement.
has a properly positioned artery mark at the bladder midpoint	Ensures proper compression of the brachial artery. Most low cost imports fail to meet this requirement from ANSI SP 10
Size Guide marking system (points out when to choose larger or smaller cuffs, properly positioned artery mark)	Prevents mis-cuffing which can lead to measurement error
Leak tested Valve	Valves are tested to ensure conformance with ANSI SP 10 for dump rate and leakage.
Premium filter equipped valve and end valve	Large valve is more ergonomic, filters prevent damaging dust build up
Unit features Latex Free bulb and bladder	Eliminates risk of latex reaction
100% inspected, assembled and packaged in the USA	Ensures unmatched consistency, quality, and performance in a value priced BP instrument
Lifetime calibration, 3 year system warranty	To protect the consumer's investment

ORDERING INFORMATION					
Schein Item#	ADC Item#	Product Description	Latex Content	U/M	Cs Pk
1126069	750W-11ABKHS	PRO PLUS 750W, Adult, Black	Latex-Free	ea.	12
1126070	752M-11ABKHS	PRO PLUS 752M, Adult, Black	Latex-Free	ea.	4