



641



647

PRIVATE LABELING	
COMPONENT	PACKAGING
N/A	2 Pc. Box/Label



PRO SPRAGUE RAPPAPORT

Our exclusive design features an improved valve mechanism that virtually eliminates acoustic leakage and a chestpiece key to facilitate after market service.

The PRO Spragues feature:

- Threaded chestpiece drum of chrome plated zinc permits the use of five interchangeable chestpiece fittings - an adult and pediatric diaphragm, and three bells (1 1/4", 1", 3/4") - for complete diagnostic capability in most patients
- Latex-free tubing
- Proprietary valve mechanism eliminates acoustic leakage
- Adjustable chrome-plated brass binaurals
- Complete with an accessory pouch containing the bell chestpieces, chestpiece key, 2 spare diaphragms, one pair of mushroom eartips, and one pair of silicon eartips
- Convenient 22" tubing length in choice of 29 fashion tubing colors
- Weighs 7 oz. (Traditional) 5.7 oz. (Single Tube)
- Overall length 30"
- 5 Year Warranty
- Inspected and packaged in the U.S.A.

Feature	Benefit
Exclusive Source	Unsurpassed quality and consistency
Chestpiece features patented valve mechanism	Virtually eliminates acoustic leakage
Convertible chestpiece can accommodate any combination of two – adult diaphragm, pediatric diaphragm, adult, pediatric, and infant bell	For the ultimate in versatility. Great scope for students, or EMS
Gasket equipped chestpiece drum	Improves leakage and improves security of chestpiece fittings
Chestpiece “key”	Facilitates removal or worn or broken diaphragms reducing the need for factory service
Dual tube design (Traditional Only)	Improves acoustic separation
Two tubing clips (Traditional Only)	Reduce incidental rubbing
Extra long binaurals	Improved practitioner comfort
Spare diaphragms	Permits self service
Includes 3 sets of eartips	For user comfort
Individually inspected and packaged in the USA	For more consistent quality
5 year warranty	To protect the consumer’s investment

ORDERING INFORMATION					
Schein Item#	ADC Item#	Product Description	Latex Content	U/M	Cs Pk
9004822	641BKHS	PRO Sprague, 22", Black	Latex-Free	ea.	20
9004823	641NHS	PRO Sprague, 22", Navy	Latex-Free	ea.	20
9004824	641NOHS	PRO Sprague, 22", Neon Orange	Latex-Free	ea.	20
9004825	641RBHS	PRO Sprague, 22", Royal Blue	Latex-Free	ea.	20
9004826	647BKHS	PRO 1-Tube Sprague, 22", Black	Latex-Free	ea.	20