

sempermed® syntegra cr

Chloroprene, Latex-Free, Fully Textured Surgical Glove

CHLOROPRENE SURGICAL Fully Textured



FEATURES

- Broad, anatomical design
- Great comfort
- Manufactured without natural rubber latex
- Fully textured
- Sterile

BENEFITS

- Reduced hand fatigue and exceptional fit
- Extraordinary latex-like fit and feel
- Eliminates exposure to latex
- Enhanced gripping
- Increased protection for use in aseptic environments

REORDER INFORMATION		Glove Color: Cream	
SCR600	6	40 PAIRS/BX	6 BX/CS
SCR650	6.5	40 PAIRS/BX	6 BX/CS
SCR700	7	40 PAIRS/BX	6 BX/CS
SCR750	7.5	40 PAIRS/BX	6 BX/CS
SCR800	8	40 PAIRS/BX	6 BX/CS
SCR850	8.5	40 PAIRS/BX	6 BX/CS
SCR-SPK		3 Pairs; Custom Sizing	

RECOMMENDED USES



SURGICAL



STERILE



OPERATING ROOM

PRODUCT SPECIFICATIONS



LATEX



TEXTURED



BEADED CUFF



SINGLE USE ONLY



HAND SPECIFIC



EXTENDED CUFF



POLYMER COATED

Our Gloves are Disposable. The People they Protect are Not.®