

SOL-CARE™ SAFETY BLOOD COLLECTION NEEDLE WITH PREATTACHED HOLDER

This safety-engineered device is intended for venous blood collection. The safety mechanism can be activated immediately after blood draw, helping to protect against needlestick injury.

#110201030006, 23 ga x 3/4" x 12", Blue
(142-3841).....25/box

#110201030002, 21 ga x 3/4" x 12", Green
(142-3842).....25/box



SOL-CARE™ SAFETY MULTI-SAMPLE BLOOD COLLECTION NEEDLE

Designed with a safety shield that completely covers and locks over the needle after use to protect against needlestick injuries.

- Safety shield aligns with a bevel up needle orientation, allowing for low-angle penetration
- Pivoting safety shield with audible click indicates it is locked in place
- Safety mechanism is designed for one-handed activation
- Safety shield is color coded to reduce potential errors

#SMSN22125, 22 ga x 1 1/4"
(142-3840) 48/box

SOL-CARE™ SAFETY BLOOD COLLECTION NEEDLE WITH MULTI-SAMPLE LUER

This safety-engineered device has a sliding safety shield that clicks when activated, alerting that the needle is safely locked.

- Transparent design allows visualization of the blood flashback and flow
- Single-handed activation option provides more flexibility to care for the patient and venipuncture site
- Sliding safety shield pushes over the needle and locks with an audible click, protecting against needlestick injury

#1102010300035, 25 ga x 3/4" x 12", Orange
(142-3845)..... 50/box

#1102010300033, 23 ga x 3/4" x 12", Blue
(142-3844)..... 50/box

#1102010300031, 21 ga x 3/4" x 7", Green
(142-3843)..... 50/box



www.henryschein.com/medical