

VK 102 Air Purifier by Radic8

Designed for Quarantine Rooms & Critical Care Environments

NOW
AVAILABLE
through
Henry Schein

Radic8 Viruskilling™ Radic8 Viruskilling™ technology was funded in 2004 by the Korean Government in response to SARS. The VirusKiller™ technology which combines air purification and photocatalytic oxidation is now in over 80% of all hospitals and schools in South Korea with over 370,000 installations.



Radic8 technology uses a two-stage process: filtration followed by photocatalytic oxidation technology.

The VirusKiller™ Technology has a **99.9999% kill rate—in a single pass**—against viruses similar to SARS-CoV-2 or harder-to-kill, such as Coronavirus DF2, Adenovirus, Influenza virus, and Poliovirus.



We Share
Clean Air

The **VK 102** patented technology combines a pre-filter, HEPA filter, activated carbon, and the reactor cell of 16 super UVC lamps with 140TiO₂ hexagon filters for coverage of spaces up to 1,776 square feet.*

VK 102 is a freestanding tower unit designed for quarantine rooms and other critical care environments. It features multi-directional airflow and can be used for negative pressure sterilization.

VirusKiller™ technology has 14 patents granted in South Korea.

Top Innovation Company Award 2018 for Infection Control as determined by InfectionControl.tips

Available Now from Henry Schein

*Assuming 8' ceiling

This product is not FDA cleared or approved. Marketed under FDA's Enforcement Policy for Air Purifiers, during the COVID-19 public health emergency. Tested on viruses that are representative of SARS-CoV-2, but not SARS-CoV-2 itself. This product has not been tested on SARS-CoV-2. This product is for medical purposes and should be used as an adjunct to personal protective equipment (PPE) and other infection control procedures. When using VirusKiller™, basic safety precautions should be observed, and instructions should be read to reduce the risk of fire, electric shock, or injury; the appliance should not be used by children or adults with reduced capabilities. Looking directly at an ultra-violet light may cause permanent damage to the eyes. Please visit <https://henryscheindigital.com/radic8>

© 2020 Henry Schein, Inc. All rights reserved. No copying without permission. Not responsible for typographical errors. 20BS9112