

# VK 103 Air Purifier by Radic8

## Designed for Open Plan Clinics, and Large Waiting Rooms

NOW  
AVAILABLE  
through  
Henry Schein

**Radic8 Viruskilling™** Radic8 Viruskilling™ technology was funded in 2004 by the Korean Government in response to SARS. The Viruskilling™ technology which combines air purification and photocatalytic oxidation is now in over 80% of all hospitals and schools in South Korea with over 370,000 installations.



**Radic8** technology uses a two-stage process: filtration followed by photocatalytic oxidation technology.

The Viruskilling™ Technology has a **99.9999% kill rate—in a single pass**—against viruses similar to SARS-CoV-2 or harder-to-kill, such as Coronavirus DF2, Adenovirus, Influenza virus, and Poliovirus.



We Share  
Clean Air

The **VK 103** patented technology combines a pre-filter, HEPA filter, activated carbon, and the reactor cell of 8 super UVC lamps with 70TiO<sub>2</sub> hexagon filters for coverage of spaces up to 1,076 square feet.\*

**VK 103** is a freestanding tower unit designed for open plan clinics and waiting rooms of all sizes, especially large rooms. It features multi-directional airflow and double filtration trays.

**Viruskilling™** technology has 14 patents granted in South Korea.

**Top Innovation Company Award 2018 for Infection Control as determined by InfectionControl.tips**

Available Now from Henry Schein

\*Assuming 8' ceiling

This product is not FDA cleared or approved. Marketed under FDA's Enforcement Policy for Air Purifiers, during the COVID-19 public health emergency. Tested on viruses that are representative of SARS-CoV-2, but not SARS-CoV-2 itself. This product has not been tested on SARS-CoV-2. This product is for medical purposes and should be used as an adjunct to personal protective equipment (PPE) and other infection control procedures. When using Viruskilling™, basic safety precautions should be observed, and instructions should be read to reduce the risk of fire, electric shock, or injury; the appliance should not be used by children or adults with reduced capabilities. Looking directly at an ultra-violet light may cause permanent damage to the eyes. Please visit <https://henryscheindigital.com/radic8>

© 2020 Henry Schein, Inc. All rights reserved. No copying without permission. Not responsible for typographical errors. 20BS9112