


# VK 401 Air Purifier by Radic8

Designed for Exam, Treatment, and Small to Medium Waiting Rooms

**NOW AVAILABLE**  
through  
Henry Schein



**Radic8 Viruskiller™ technology** was funded in 2004 by the Korean Government in response to SARS. The VirusKiller™ technology which combines air purification and photocatalytic oxidation is now in over 80% of all hospitals and schools in South Korea with over 370,000 installations.

Radic8 technology uses a two-stage process: filtration followed by photocatalytic oxidation technology. The Viruskiller™ Technology has a **99.9999% kill rate—in a single pass—against viruses** similar to SARS-CoV-2 or harder-to-kill, such as Coronavirus DF2, Adenovirus, Influenza virus, and Poliovirus.

The **VK 401** patented technology combines pre-filter, HEPA filter, activated carbon, and the reactor cell of 8 super UVC lamps with 40TiO<sub>2</sub> hexagon filters for coverage of spaces up to 645 square feet.\*

With its large air exchange capacity, the **VK 401** is designed for exam rooms, treatment rooms, and small to medium waiting rooms. It can be wall-mounted to save space or freestanding. A convenient remote control and timer are included.

**Viruskiller™** technology has 14 patents granted in South Korea.

**Top Innovation Company Award 2018 for Infection Control as determined by InfectionControl.tips**

Available Now from Henry Schein



\*Assuming 8' ceiling.  
This product is not FDA cleared or approved. Marketed under FDA's Enforcement Policy for Air Purifiers, during the COVID-19 public health emergency. Tested on viruses that are representative of SARS-CoV-2, but not SARS-CoV-2 itself. This product has not been tested on SARS-CoV-2. This product is for medical purposes and should be used as an adjunct to personal protective equipment (PPE) and other infection control procedures. When using Viruskiller™, basic safety precautions should be observed, and instructions should be read to reduce the risk of fire, electric shock, or injury; the appliance should not be used by children or adults with reduced capabilities. Looking directly at an ultra-violet light may cause permanent damage to the eyes. Please visit <https://henryscheindigital.com/radic8>  
© 2020 Henry Schein, Inc. All rights reserved. No copying without permission. Not responsible for typographical errors. 20BS8610