

# Z-800F

## IV Infusion Pump

**ZYNO SOLUTIONS**  
Intelligent Infusion Systems



*Simple*  
*Smart*  
*Safe*  
*Dependable*

### Infusion Pump System

- Infusion technology for oncology, long term care, specialty pharmacy and other disciplines
- Designed with direct input from nurses
- Robust performance and reliability
- Reduced repair issues creates less downtime
- Upgradeable software for continuing feature enhancements

# Z-800F IV Infusion Pump

Simple. Safe. Smart. Dependable

## Simplified Set-up

- Simple 3 step intuitive programming
- Easy to load straight line IV sets
- Reliable long lasting quick charge battery

## Reliable Clinical Solutions

- Multi-therapy infusion mode options
- Extended nurse light notification
- Set based anti-free-flow protection

## Customized Infusion Workflow

- Quick store and retrieve protocols
- Simple customizable mode options
- 10 step infusion for multi-step infusions



## Specifications

Mechanism	Linear Peristaltic	Power	100-240 VAC~, 50-60Hz
Flow Rate Range	1-1200 mL/Hr	Leakage Current	≤100uA
Volume To Be Infused (VTBI)	1-9,999 mL	Battery	9.6V, 4.7Ah NiMH (rechargeable)
Accuracy	+/- 5%	Battery Operating Time	8 Hours at 125 mL/hr
Infusion Modes	Continuous Rate/Volume, Continuous Time/Volume, Blood, Ten Step Sequencing, TPN, Protocol, Intermittent	Air-in-Line	Ultrasonic, 14 level Adjustable
Free Flow Protection	Pump based protection on all models; added IV set based protection in Z-800F and Z-800WF	Occlusion Pressure	4-30 PSI Adjustable
Housing	Metal	Tubing	Zyno standard tubing is not made with natural rubber latex, and is not made with DEHP
Weight	7.0 lbs (3.2 kg)	Alarms	Air-in-Line, Occlusion, Door Open, Battery Empty, System Error
Size	8.6 x 5.7 x 5.3 in (HxWxD)	Warnings	Infusion Complete - KVO, Pump Unattended, Near End, Low Battery

www.zynosolutions.com | info@zynosolutions.com | 177 Pine Street, Natick, MA 01760