

NEW YEAR, NEW PRODUCTS, NEW DEALS!



SensiLase® PAD-IQ™


orthopedics
equipment
pharmaceuticals



20/20 Imaging



HENRY SCHEIN®
BRAND PRODUCTS



HENRY SCHEIN®
BRAND PRODUCTS



Alma Lasers™
Wellbeing Through Technology

MENTION
PROMO
CODE **4F** TO ENSURE
PROMOTIONAL
PRICING

infection control

instruments

FOOT & ANKLE



JOIN OUR COMMUNITY

To Order: 1-800-323-5110 8am-9pm (et)

To Fax: 1-800-524-4969 24 Hours

www.henryschein.com/footandankle

Celebrating 80 Years of Excellence



Since Henry Schein opened his pharmacy in Queens, New York in 1932, our organization has become a Fortune 300 company and the largest provider of health care products and services to office-based practitioners. With more than 14,000 Team Schein Members and operations or affiliates in 25 countries, Henry Schein, Inc. remains committed to serving its customers, supplier partners, investors, society, and Team Schein.

Your Success Is Our Goal

Resale products, equipment, and the utilization of practice solutions and EHR technology are key components for generating increased revenue. Henry Schein Medical provides the Foot and Ankle Specialist with medical supplies, equipment, pharmaceuticals, and technology solutions that will support business goals and enhance patient services. Whatever you need can be found in our latest sale flyer or online, along with the personal service your practice demands.

Visit our Web site and take advantage of features like inventory management and billing history: www.henryschein.com/footandankle.



PRIVILEGES Loyalty Rewards Program

Take advantage of our PRIVILEGES Loyalty Rewards Program and earn points for all of your online orders. Plus, you can earn up to 115,000 bonus points with your Value Certificates. Sign up today—enrollment is fast and easy!

For more information, visit the Web site above and click on the PRIVILEGES bar/banner on the home page or visit: www.henryschein.com/privilegesmd.

Henry Schein Brand Products



Look for the Henry Schein Brand products throughout this flyer. You can rely on the Henry Schein Brand name for high-quality products that provide extra value at reasonable prices. Additionally, our product line has met our high standards to earn the Henry Schein Seal of Excellence—**your guarantee of satisfaction.**

Henry Schein Medical is your reliable network of trusted advisors.

Rely on Us for...

Expert Guidance

Our knowledgeable team of trusted advisors is dedicated to collaborating with you to meet specific goals related to quality patient care, staff and team growth, and profitability.

- Revenue Enhancement Tools
- Cost-Containment Programs
- Staff Training and Team Building
- Practice Marketing
- Technology Solutions

Price Integrity

We have relationships at the local, regional, and national level—including all major GPOs—to help manage operating costs and drive value for your practice. Web-enabled reporting tools guarantee your pricing is accurate and consistent.

Guaranteed Fast Delivery

All orders received by 5 PM local time are shipped the same day. Your products arrive within 24–48 hours.

Pharmaceuticals and Diagnostics

Our pharmaceutical and diagnostic solutions include brand names as well as a complete line of generics and Henry Schein Brand options for added value.

TABLE OF CONTENTS

equipment	4–9
orthopedics	10–14
instruments	15–18
infection control ...	19–20
pharmaceuticals ...	21–23

 **HENRY SCHEIN**[®]
MEDICAL | FOOT & ANKLE

Rely on Us

Alma® Lasers
Now Offers the
Only Platform with
**Two Distinct Laser
Technologies** for
Onychomycosis*
(nail fungus)

- Fast, versatile, easy to use
- Visible results
- Minimal discomfort, no downtime, topicals or disposables

Discover ClearChoice and how it can improve patient satisfaction and help increase revenue for your practice.



 **ClearChoice™**
Complete Nail Restoration



HARMONY® ClearChoice™

Contact your Henry Schein Representative to order

Alma® Lasers
Wellbeing Through Technology™

www.almalasers.com

HENRY SCHEIN®
MEDICAL | FOOT & ANKLE

* FDA cleared for use for the temporary increase of clear nail in patients with onychomycosis.

mention
promo
code **4F** to ensure
promotional
pricing

www.henryschein.com/footandankle

13MS7514



Peripheral Vascular Testing Just Got Smarter!

Introducing...

SensiLase® PAD-IQ™



- Complete Non-Invasive Lower Extremity Vascular Assessment
 - Skin Perfusion Pressure - Microcirculation (Reactive Hyperemia)
 - Pulse Volume Recording
 - Ankle Brachial Index (ABI) Calculation
 - Peripheral pressure measurement
- Simultaneous Bilateral Testing
- Touch screen, no mouse or keyboard needed. Easy to use.
- Fully automated operation
- Portable: Table, IV pole or bed mount
- Full communications and connectivity for customizable data import/export/documentation

SensiLase® PAD-IQ™ The Intelligent Vascular Test

Exclusive to Henry Schein Podiatrists, item # 614-3004



When you want the best for your practice, We Have Your Solution at MTI.

MTI 527S/P SERIES PODIATRY CHAIR FEATURES:

POWER LIFT, BACK, AND TILT FUNCTIONS

Superior positioning for all your patient needs

SLIDE BACK REMOVABLE FLOATING ARMS™

Manage patient arm positioning with comfort and ease

SUPERIOR LIFT CAPACITY

Lifts up to 650 lbs effortlessly

NEW! REMOVABLE LEG WRAPPING SUPPORT*



Foot/Heel wrapping position



Leg wrapping position

For more information and incentives on our products:

Call (800) 924-4655

Email sales@mti.net

Visit our web site at www.mti.net

*Patent pending

Strength in patient care.™



Combined weight of Brad, Jeff, and Mike is 505lbs.

3655 Ninigret Drive, Salt Lake City, Utah 84104
800.924.4655 • 801.875.4999 • Fax 801.952.0548

www.mti.net



6



HENRY SCHEIN® PRO PLUS SERIES CLOCKFACE ANEROIDS

Each Aneroid features:

- Large, easy-to-read 6" luminescent dial
- 300mm-Hg, CE marked, no-pin stop movement
- Inflation system includes black nylon cuff, latex inflation bladder, bulb and valve, and 8' length

Wall Aneroid includes:

- Wall bracket that swivels through 110° rotation for easy viewing with integral cuff storage

Mobile Aneroid features:

- Height-adjustable Spider-Leg™ caster stand with range of 34" to 48" (from bottom of base to top of unit)
- Integral ABS cuff-storage compartment
- Manometer is angled back 15° for easy viewing
- 18" base diameter

Wall Aneroid
Latex-Free
(112-6069)ea

Mobile Aneroid
Latex-Free
(112-6070)ea

WelchAllyn

ProBP™ 3400

Small enough to carry everywhere, with the same accuracy as our vital signs devices, the Connex® ProBP™ 3400 digital blood pressure device is a smart choice for clinical settings. Take fast, accurate blood pressure readings in 15 seconds, while USB-wired connectivity lets you send them straight to your EMR/HIS. Includes medium cuff. Not intended for neonatal use.

With SureBP™ Thermometry
Handheld, #34XFHT-B
(566-0229)ea

With Standard BP
Handheld, #34XXHT-B
(566-0227)ea

With SureBP™ Thermometry
Wall Mount, #34XFWT-B
(566-0233)ea

With Standard BP
Wall Mount, #34XXWT-B
(566-0232)ea

With SureBP™ Thermometry
Mobile Stand, #34XFST-B
(566-0238)ea

With Standard BP
Mobile Stand, #34XXST-B
(566-0237)ea



Ask Us About Our Great Leasing Programs!



iCR3600 SF Foot X-ray Package

- High-resolution images from phosphor plate technology is coupled with adaptive image processing software
 - Provides customized podiatry labels for faster searching of images and communication, specific podiatric measurements such as IMA and HAA, and customizable templates for preop/postop reporting and communication
 - Cassette sizes: 10" x 12" (25 cm x 30 cm), 8" x 10" (20 cm x 25 cm)
 - Image access time: 25 seconds (10" x 12")
 - Plate throughput: 83 plates per hour (10" x 12")
 - Dimensions: 18"H x 34.5"W x 10.5"D (45.7 cm H x 87.6 cm W x 25.4 cm D)
 - Wall-mountable or tabletop
- (116-2133).....ea**



20/20 Digital Imaging System

FEATURES AND BENEFITS:

- Patented Tilt Head System (eliminates need to rotate patient for AP views)
- 3½" high base (improved height for patient use)
- Patient support rail allows access from side or front
- 29" SID AP/oblique view
- 36" SID lateral view
- "SR" series proprietary tube
- Low-voltage exposure switch
- Qualifies for ADA tax credit

SPECTACLE SOFTWARE:

- Innovative "Point, Click & View" process
 - Cleaner artifact-free data
 - Advanced DICOM image processing
 - 20/20 Image Management Software
 - iSPECTACLE web-based software (optional) includes:
 - Installation and in-service applications
 - Full 12-month parts and labor warranty
 - 24/7 technical support
- (115-6770)ea**



Pro Vac-File Drill

Powerful yet quiet vacuum is built into the ergonomic handpiece. With variable and reversible speed up to 40,000 rpm, the Pro Vac File makes quick work of most debriding tasks. Lightweight (only 5 oz) handpiece fits comfortably in your hand and has its own On/Off switch. The auto-locking chuck eliminates the need for wrenches or twisting when changing tools. Filtration is .02 microns at 99%. 2-year warranty. Foot pedal included.

(115-6992)ea



Nail Glide 2100 Super Drill

- 35,000 rpm
- High torque
- Lightweight (6 oz), autoclavable handpiece
- No noise
- No vibration—smooth operation
- Forward and reverse
- Locking chuck
- Variable-speed foot pedal
- 1-year warranty with an optional 2-year warranty

(158-0002).....ea



Pro Power Debriding Drill

- 30,000 rpm
- High torque
- Lightweight (4 oz), autoclavable handpiece
- Lightweight power station
- 10-hour battery life
- No noise
- No vibration—smooth operation
- Forward and reverse
- Locking chuck
- 1-year warranty with an optional 2-year warranty

(158-0003).....ea

Functionality, efficiency and versatility - Midmark 647.

With its best-in-class foot section design, the Midmark 647 is the next step forward in foot, ankle and wound care. Innovative and ergonomically advanced, the 647 helps you work at peak performance, even during lengthy procedures.

- Integrated touch pad controls for convenient point of care access
- Ergonomic, conveniently positioned release handles provide smooth, one-handed operation of the foot section
- Safe for all patients. Designed for podiatry demographics, the 647 has the size and stability to safely position patients for podiatry examination and treatments.

For more information, call **1.800.323.5110** or visit **midmark.com/ScheinQ1FootAnkle**.



Midmark Corporation, Versailles OH USA



Ritter M11 UltraClave® Automatic Sterilizer

Like the Ritter M9, this unit also contains an automatic opening door. Unique chamber size gives it the most sterilizing capacity of any standard countertop sterilizer on the market. Featuring a modern look, quiet door operation, and pulse air-removal technology, this sterilizer couldn't be easier, more efficient, or more reliable. You load it, set it, and forget it.

- Door opens automatically
- UL listed and CSA/ASME certified
- Chamber size: 11" x 18"
- Weight: 99 lb

(386-8917)ea



M9 UltraClave & Sterilizer

This reliable and efficient sterilizer features a sleek look, quiet-door operation, and pulse-air-removal technology. Front fill-and-drain tube feature means you never have to move it to maintain or drain the water. Load it, set it, and forget it.

- Reliable and easy to use
- Provides quality infection control
- Select your cycle and time is set
- Press "Start" and M9 fills with correct amount of water
- When cycle is complete, door opens automatically and quietly for fast drying time

(386-7984)ea

Histofreezer[®]

Portable Cryosurgical System

A Cool Way to Remove Common Warts & Benign Lesions



CLINICALLY PROVEN, SAFE AND EASY TO USE

- Results are comparable to Liquid Nitrogen
- Patient friendly, low risk of pain and scarring
- No Freeze-Thaw-Freeze procedure
- Economically priced and reimbursable
- 3-year shelf-life

<p>Histofreezer 60M2C (810-0004) Mfg. #1001-0292</p> <ul style="list-style-type: none"> • 60-120 applications per cryogen kit • (36) 5mm + (24) 2mm applicators 	<p>Histofreezer 36M1C (810-0005) Mfg. #1001-0295</p> <ul style="list-style-type: none"> • 36-72 applications per cryogen kit • (24) 5mm + (12) 2mm applicators
<p>Histofreezer 40T1C (810-0006) Mfg. #1001-0293</p> <ul style="list-style-type: none"> • 40-80 applications per cryogen kit • (40) 2mm applicators 	<p>Histofreezer 36F1C (810-0007) Mfg. #1001-0294</p> <ul style="list-style-type: none"> • Ideal for Podiatry and larger lesions • 36-72 applications per cryogen kit • (36) 5mm applicators



1. Data on file, OraSure Technologies, Inc., Bethlehem, PA.
© 2012 OraSure Technologies, Inc. All Rights Reserved. HIST00331 (rev. 10/12)

Go to www.Histofreezer.com/us/ to order your **FREE** sample kit today!



DARCO

Body Armor® Vario

This brace is a general-use ankle brace that adds an element of customization. The brace is designed for use in cases of PTTD, tendonitis, lateral and medial ankle sprains, and ligament issues, as well as a corrective and postoperative device. This brace comes with 3 posting discs of 2, 4 and 6 degrees to correct varus/valgus alignment, and 3 arch supports—high, medium and low—to provide the support needed allowing for customized application. This brace is comfortable and comes with a neoprene sleeve for additional compression and support if needed. Black.

- Ankle Braceea
- Specify:
- Left(765-0008)
- Right.....(765-0009)



DARCO

DARCO FX Pro™ Low and High Walkers

Low-profile walkers that match the heel height of most shoes. They are very light, weighing less than 2 lb. each.

- Low**
- X-small (240-9624)ea
- Small (240-9497)ea
- Medium (240-4960)ea
- Large (240-4658)ea
- X-large (240-5795)ea
- High**
- X-small (240-3691)ea
- Small (240-3642)ea
- Medium (240-2969)ea
- Large (240-0084)ea
- X-large (240-2946)ea



DARCO

Body Armor Stirrup Ankle Brace

Lightweight, low-profile, one-size-fits-all ankle product that provides excellent compression and motion restriction. The anatomically contoured stirrups and soft neoprene lining deliver a combination of rigidity and comfort. This is a true convertible brace, fitting both left and right ankles.

- Universal Size (116-9812).....ea



DARCO

DARCO All-Purpose Boots

- Black
- X-Small (948-3668) #APQ0B.....ea
- Bootsea
- Specify:
- #APQ1B, Small(104-3408)
- #APQ2B, Medium(104-3409)
- #APQ3B, Large(869-6446)
- #APQ4B, X-large.....(941-3537)



DARCO

DARCO Softie Shoes

- Men's Small (812-6086)ea
- Men's Medium (812-4388)ea
- Men's Large (812-3658)ea
- Women's Small (812-6517)ea
- Women's Medium (812-8848)ea
- Women's Large (812-6799)ea



DARCO

Body Armor® Dorsal Night Splint

- One size fits all. (113-7635).....ea



DARCO

DARCO Pneumatic Night Splints

- Nonskid sole.
- Small (Up to: Men's 4½, Women's 6½) (685-6907)ea
- Medium (Men's 5-8, Women's 7-10) (685-2884)ea
- Large (Men's 8½-11, Women's 10½-13) (685-7200)ea
- X-large (Men's 11½-14, Women's 13½+) (685-0383)ea



SWED-0, INC.
Strap Lok®

Ultrathin and Lightweight
Fits right or left ankle.

- Braceea
- Specify:
- X-small, White(689-5199)
- Small, White(364-3028)
- Medium, White(364-3204)
- Large, White(364-9493)
- X-large, White(689-8387)
- XX-large, White(487-4150)
- X-small, Black(101-1497)
- Small, Black(364-4040)
- Medium, Black.....(364-7196)
- Large, Black(364-4225)
- X-large, Black(364-0817)
- XX-large, Black.....(674-2199)



ProTech Control Full-Length Orthotics

Unexcelled relief from heel and arch pain, as well as incredible support, control, and stability. Features Heel Secure System to enhance foot control, 2" extrinsic post that stabilizes the foot, deeper heel cup to better protect the heel, double strength arch firmness, double layer of Poron & EVA cushioning, Poron heel pad for added comfort, antimicrobial top cover for dryness. Maximum support and maximum comfort. Fits casual and athletic shoes, motion control and stability running shoes. Available in 11 sizes.

Orthoticspair

Specify:

- #PCTa, Men's 3-3½; Women's 5-5½(118-0098)
- #PCTA, Men's 4-4½; Women's 6-6½(118-0099)
- #PCTB, Men's 5-5½; Women's 7-7½(118-0100)
- #PCTC, Men's 6-6½; Women's 8-8½(118-0101)
- #PCTD, Men's 7-7½; Women's 9-9½(118-0102)
- #PCTE, Men's 8-8½; Women's 10-10½(118-0103)
- #PCTF, Men's 9-9½; Women's 11-11½(118-0104)
- #PCTG, Men's 10-10½; Women's 12-12½(118-0105)
- #PCTH, Men's 11-11½; Women's 13-13½(118-0106)
- #PCTIJ, Men's 12-13½; Women's 14-15½(118-0107)
- #PCTK, Men's 14-14½(118-0108)



Powerstep Protech ¾ Orthotics

¾ length prescription-like arch support calibrated with flexibility. Heel cradle with opening filled with Poron cushioning, resilient cushioned casing, antimicrobial top cover. Fits a broad range of shoes. Available in 7 sizes.

Protech ¾ Orthoticspair

Specify:

- #P2ASM, Men's 4-4½, Women's 6-6½(108-1324)
- #P2A, Men's 5-6½, Women's 7-8½(108-0298)
- #P2B, Men's 7-8½, Women's 9-10½(108-1328)
- #P2C, Men's 9-9½, Women's 11-11½(108-1329)
- #P2D, Men's 10-10½, Women's 12-12½(108-1330)
- #P2E, Men's 11-11½, Women's 13-13½(108-1332)
- #P2F, Men's 12-13½, Women's 14-15½(108-1334)



littleSTEPS™ Foot Orthotics

- UCBL-type design with 30-mm heel cup depth to maximize rear foot control and realign the Achilles tendon
- Designed with medial and lateral flanges, medial rear foot posting, and skive
- Arch support with a deep flare to ensure both comfort and correction
- A unique thermoplastic compound, which has a combination of strength and a "soft-edge" feel for patient comfort
- Durable and easy to clean with soap and water
- Fits in most standard athletic shoes
- Treats Sever's disease, flexible pes planus, and idiopathic toe-in or toe-walking gait

Foot Orthoticspair

Specify:

- Size 1, Toddler's 8½-9½(125-0117)
- Size 2, Toddler's 10-11(125-0119)
- Size 3, Toddler's 11½-12½(125-0120)
- Size 4, Toddler's 13-Children's 1(125-0121)
- Size 5, Children's 1½-2½(125-0122)
- Size 6, Children's 3-4(125-0123)
- Size 7, Children's 4½-6(125-0124)



Toe Crest with Reinforced Toe Loop



Prevents callosities and hammer toes. Stretches contracted toes. Helps reduce pressure and pain caused by hammer, claw, and mallet toes.

Toe Crestsea

Specify:

- Right, Small(114-8210)
- Left, Small(114-8215)
- Right, Medium(114-8218)
- Left, Medium(114-8219)
- Right, Large(114-8221)
- Left, Large(114-8222)



Gel Digital Pad Medpack



Digital pad on a tube. Soothes lesions caused by pressure or rubbing and reduces pain caused by calluses, blisters, soft corns, and hammer toes.

Small/Medium (114-1136)12/pkg

Large/X-large (114-1134)12/pkg



Gel Hallux Straightener



Helps realign hallux valgus deformities by preventing the big toe from overlapping onto the second. Prevents and reduces pain caused by calluses and ingrown toenails.

(114-4085) Small4/pkg

(114-4089) Medium4/pkg

(114-4091) Large4/pkg



GelSmart Mesh Digital Caps

Mesh fabric fully coated with soft gel. Ideal for hammer toes, corns, calluses, nails, and cuticles.

Small/Medium (108-2819)6/pkg

Large/X-large (112-3537)6/pkg



Achilles Gel Heel Protection Sleeve

Provides gentle compression. Ideal for pump bumps; Haglunds deformity; and sore, fatigued or injured Achilles tendon.

Small/Medium (108-2617)ea

Large/X-large (108-2619)ea



GelSmart Ribbed Digital Caps

Flexible fabric half-coated with soft gel. Ideal for corns, blisters, and other injuries to digits.

Small/Medium (110-5622)6/pkg

Large/X-large (110-5614)6/pkg



Arch Binder with Metatarsal Pad

Supports and aligns metatarsal heads.

Small (634-0513)pair

Medium (102-1295)pair

Large (797-5852)pair



Turf Toe Straps

Soft, heavy-duty moleskin with adhesive backing. Ideal for turf toe, hyperextended toe, and pre- and post-abduction of the hallux.

#TT61, 7⁵/₈" x 4"
(399-4556)100/bag



Tubular Foam Sleeves

Soft, breathable foam lined with seamless fabric.

- Reduces friction and pressure
 - 12" long; cut to desired length
- 1" Overlap
(335-6790)8/pkg
- 5/8" Single
(843-1657)12/pkg



Orthopedic Felt (Adhesive)

- 6" x 2¹/₂-yd Rolls
- 1/8" Thick
(242-9063)ea
- 1/4" Thick
(242-8317)ea



Felt, Foam, and Moleskin Rolls

Made of high-quality orthopedic felt, foam, and moleskin, these materials are used to customize padding for injured area, in the manufacture of splints, and in blister care. Adhesive backed.

- Moleskin, 9" x 4 yd
(155-2385)ea
- Moleskin, 12" x 5 yd
(379-6681)ea
- Moleskin, 22" x 5 yd
(242-5879)6 rolls/box
- 1/4" Thick White Foam, 6" x 2¹/₂ yd
(242-2969)ea
- 1/8" White Felt, 6" x 2¹/₂ yd
(155-1527)ea
- 1/4" White Felt, 6" x 2¹/₂ yd
(155-2461)ea
- 1/8" Flesh Felt, 6" x 10 yd
(612-3765)ea



Pedifoam Single

- Single
3/4" Medium
(962-3888)10/Pk
- 1"
Large
(335-3573)8/PK



SEAL-TIGHT® Original

Cast and Bandage Protectors, Adult Size
Recommended by doctors for therapeutic care and daily showering and bathing. Easy to apply with the patented application ring, which eliminates the need for strapping, pumps, or Velcro® closure. Nonlatex polyvinyl bag is durable, reusable, and stretches easily over cast or bandage to form a watertight seal. With normal daily care, will last 6-8 weeks.

- Legs
Foot/Ankle, 12" (28 cm)
(637-2389)ea
- Short Leg, 24" (59 cm)
(802-0475)ea
- Wide Short Leg, 23" (58 cm)
(838-0829)ea
- Long Leg, 43" (106 cm)
(948-0453)ea

Henry Schein has all your casting products covered



HENRY SCHEIN®

FIBERGLASS CASTING TAPE, 4-yd Rolls

- | | |
|------------------------|------------------------|
| 2"10 rolls/box | 3"10 rolls/box |
| Specify: | Specify: |
| White(900-4635) | White(900-4636) |
| Blue(900-4638) | Blue(900-4639) |
| Black(900-4641) | Black(900-4642) |
| Pink(900-4644) | Pink(900-4645) |
| Red(900-4646) | Red(900-4647) |
| Purple(900-4648) | Purple(900-4649) |
| Green(900-4650) | Green(900-4651) |

- 4"10 rolls/box
- Specify:
- White(900-4637)
- Blue(900-4640)
- Black(900-4643)



HENRY SCHEIN®

SYNTHETIC CAST PADDING

- Non-Sterile,
4-yd Rolls
- 2"
(900-4997) ...24 rolls/bag
- 3"
(900-4998) ...12 rolls/bag
- 4"
(900-4999) ...12 rolls/bag
- 6"
(900-7001)6 rolls/bag



HENRY SCHEIN®

SYNTHETIC STOCKINETTES

- Non-Sterile, 25-yd Rolls
- 2"
(112-5513)ea
- 3"
(112-5514)ea
- 4"
(112-5515)ea
- 6"
(112-5516)ea



HENRY SCHEIN®

COTTON CAST PADDING

- Non-Sterile, 4-yd Rolls
- 2"
(900-7023)24/bag
- 3"
(900-7024)12/Bag
- 4"
(900-7025)12/bag
- 6"
(900-7026)6/bag



PROCARE

STEP UP TO VALUE



An affordable and comprehensive family of Walking Braces, the **XCELTRAX™** provides ideal support and comfort to help promote healing for patients with mild to severe foot and ankle injuries.

HSI Item Code	Product Description
1530240	XcelTrax Air, Ped/XS
1530241	XcelTrax Air, Small
1530242	XcelTrax Air, Medium
1530243	XcelTrax Air, Large
1530244	XcelTrax Air, X-Large
1530246	XcelTrax Air Ankle, Ped/XS
1530247	XcelTrax Air Ankle, Small
1530248	XcelTrax Air Ankle, Medium
1530249	XcelTrax Air Ankle, Large
1530250	XcelTrax Air Ankle, X-Large

HSI Item Code	Product Description
1530251	XcelTrax, Ped/XS
1530252	XcelTrax, Small
1530253	XcelTrax, Medium
1530254	XcelTrax, Large
1530255	XcelTrax, X-Large
1530258	XcelTrax Ankle, Ped/XS
1530259	XcelTrax Ankle, Small
1530262	XcelTrax Ankle, Medium
1530263	XcelTrax Ankle, Large
1530266	XcelTrax Ankle, X-Large

For more information on the entire product offering of XcelTrax, contact your Henry Schein Representative. Or to order call: 1-800-772-4346; fax: 1-800-329-9109; or visit us at www.henryschein.com/medical



DeRoyal® Element® Sport

Powered by The Boa® Closure System

For more information, application videos and product detail, visit www.derooyal.com/medicalproducts/boa.aspx.



MICRO-ADJUSTABILITY | GLOVE-LIKE FIT | ON-THE-FLY ADJUSTMENT



Element® Sport - L1971		
SIZE	RIGHT	LEFT
Small	3720095	3720098
Medium	3720096	3720099
Large	3720097	3720100

Shoe Size Information		
SIZE	MEN	WOMEN
Small	< 8	< 9½
Medium	8 - 12	9½ - 13½
Large	12 +	13½ +

The **ELEMENT® SPORT** is for severe ankle sprains and overuse issues such as tendonitis.



DeRoyal
Improving Care. Improving Business.

NEW!

CoFlex™ • UBC Kit CoFlex™ • UBZ Kit

... a 2-step level of comfort and healing unlike any other Unna boot on the market!



Andover's Unna Boot Kit	vs.	Currently in the market
• Foam		• Gauze
• Absorbs 20x's its dry weight		• Absorbs 2.7x's its dry weight
• Zinc/Calamine evenly distributed		• Uneven distribution of Zinc/Calamine
• Clean, no mess to clean up		• Messy for clinician
• No indentation/skin breakdown		• Creates indentation due to pressure
• Stays moist		• Many dry out within the week
• Delivers consistent compression /prevents edema		• Inconsistent compression
• Comes in a kit, including a stocking		• Rolls sold separately
• Qualifies for 2 HCPCS Codes: A6456 (Step 1) & A6454 (Step 2)		• Only qualifies for 1 HCPCS Code: A6456



Calamine		Zinc	
Item #6811621 3" Unna Boot Kit Latex	Item #6811632 3" Unna Boot Kit Latex-Free	Item #6811629 3" Unna Boot Kit Latex	Item #6811635 3" Unna Boot Kit Latex-Free
Item #6811638 4" Unna Boot Kit Latex	Item #6811651 4" Unna Boot Kit Latex-Free	Item #6811643 4" Unna Boot Kit Latex	Item #6811656 4" Unna Boot Kit Latex-Free

U.S. PATENT NO. 5,762,623 | U.S. PATENT NO. 6,156,424 | EUROPEAN PATENT NO. 0839222 | PATENT PENDING

For a **FREE** sample of CoFlex UBC/UBZ, call 800.432.6686 or email sample@andoverhealthcare.com
Please reference Henry Schein Foot & Ankle Ad.



Proudly Made in USA



3M™ Coban™ LF Latex Free Self-Adherent Wrap with Hand Tear



- Now with hand tear
- Soft to the touch
- For securement, compression, immobilization & more



Clinicians like you rock!

You know it, your patients know it and we know it. That's why we've made your brand of choice self-adherent wrap even easier to use.

Now with a scissors-free application, as well as the soft feel and dependable hold you expect, Coban LF Wrap with Hand Tear lets you tear it up!

Contact your Henry Schein representative for pricing.

For securement, compression, immobilization and all the ways you use it!

Tear it up!



3M and Coban are trademarks of 3M. © 3M 2012. All Rights Reserved.



HENRY SCHEIN®
EXCAVATOR CURETTE 5½"
 German, Double End, Small 1.5mm and Medium 2mm
 (104-5536)ea

HENRY SCHEIN®
NAIL SPLITTER 5"
 English Anvil Pattern, 5" Double Spring, Stainless Steel
 (104-8507)ea

HENRY SCHEIN®
NAIL SPLITTER
 German, Straight Jaws, Double Spring, Stainless Steel
 4" (104-5512)ea
 5" (104-5526)ea
 5", Extra Narrow (104-7539)ea

HENRY SCHEIN®
DOUBLE-ACTION NAIL NIPPER 6"
 German, Sharp, Concave Jaws, Single Barrel, Stainless Steel
 (104-5530)ea

HENRY SCHEIN®
NAIL NIPPER 5½"
 German, Concave Jaws, Double Spring, Stainless Steel
 (102-5405)ea

HENRY SCHEIN®
NAIL NIPPER 5½"
 German, Concave Jaws, Barrel Spring, Chrome
 (102-6341)ea

HENRY SCHEIN®
RUSKIN BONE-CUTTING FORCEPS 6"
 German, Standard Blades
 (104-6986)ea

HENRY SCHEIN®
BROWN-ADSON DRESSING FORCEPS 4¾"
 German, 7 x 7 Side-Grasping Teeth
 (101-6911)ea



HENRY SCHEIN®
ADSON DRESSING FORCEPS 4¾"
 Serrated, Delicate Pattern
 (100-5659)ea



HENRY SCHEIN®
ADSON TISSUE FORCEPS 4¾"
 German, 1 x 2 Teeth
 (101-4803)ea



HENRY SCHEIN®
ADSON DRESSING FORCEPS 4¾"
 German, Serrated, Delicate Pattern
 (101-6279)ea



HENRY SCHEIN®
KELLY HEMOSTAT 5½"
 German Stainless Steel, Straight
 (100-3357)ea



HENRY SCHEIN®
KELLY FORCEPS 5½"
 German Stainless Steel, Curved
 (100-8381)ea



HENRY SCHEIN®
HALSTEAD MOSQUITO FORCEPS 5"
 German Stainless Steel, Straight
 (100-0846)ea



HENRY SCHEIN®
WEITLANER RETRACTOR 4"
 German, Blunt Points, 2 x 3 Teeth
 (104-9635)ea



HENRY SCHEIN®
LISTER BANDAGE SCISSORS 5½"
 German, Extra Fine
 (104-8132)ea



HENRY SCHEIN®
LISTER BANDAGE SCISSORS 5½"
 German, Extra Fine
 (100-6589)ea



HENRY SCHEIN®
LISTER BANDAGE SCISSORS 7¼"
 German, Extra Fine
 (100-3334)ea



HENRY SCHEIN®
LISTER BANDAGE SCISSORS 5½"
 Extra Fine
 (100-1087)ea



HENRY SCHEIN®
BANDAGE AND UTILITY SCISSORS 7"
 German
 (104-9280)ea



HENRY SCHEIN®
OPERATING ROOM SCISSORS 5½"
 German Stainless Steel, Straight Scissorsea
 Specify:
 Sharp Points.....(100-4549)
 Sharp/Blunt Points(100-7561)



HENRY SCHEIN®
SPENCER STITCH SCISSORS 3½"
 German Stainless Steel
 (100-6024)ea



HENRY SCHEIN®
DISPOSABLE BIOPSY PUNCH
 Available in 9 different sizes from 1.5 mm to 8.0 mm. Sterile, individually pouched with a tamper-evident seal. Our disposable skin biopsy punches have beneficial properties not found with generic products. Features include:
 • Stainless steel ground blade
 • Seamless punch
 • Ergonomic blue handle
 Biopsy Punch.....ea
 Specify:
 1.5 mm.....(900-4863)
 2 mm.....(900-4864)
 2.5 mm.....(900-4865)
 3 mm.....(900-4866)
 3.5 mm.....(900-4867)
 4 mm.....(900-4868)
 5 mm.....(900-4869)
 6 mm.....(900-4870)
 8 mm.....(900-4871)


Integra® Miltex®

Podiatric Instrumentation

Limit uncertainty with high quality instrumentation for your everyday surgical and podiatric needs

- Nippers
- Chisels
- Splitters
- Elevators
- Curettes
- Disposables

INTEGRA
LIMIT UNCERTAINTY



For more information or to place an order, please contact:
Phone: 800-P-SCHEIN (800-772-4346)
Fax: 800-329-9109
Website: www.henryschein.com/medical

Illustrations and content provide general description only and may be subject to change.
Integra, the Integra logo, and Miltex are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries.
Copyright © 2012 Integra LifeSciences Corporation
Schein Medical Foot and Ankle Ad Q4/2012 PM11120

**SUTURES
ABSORBABLE**

All Sutures..... 12/Box

ITEM CODE	NEEDLE	LENGTH	SIZE	MATERIAL
100-6120	C-9	18"	3-0	PLAIN GUT
101-3192	C-3	18"	5-0	PLAIN GUT
101-2976	C-6	27"	3-0	PLAIN GUT
101-0668	C-6	27"	4-0	PLAIN GUT
100-7771	C-6	18"	4-0	PLAIN GUT
900-7471	PC-1	18"	5-0	PLAIN GUT
900-7485	DSCM-13	18"	6-0	PLAIN GUT
101-1582	C-6	27"	3-0	CHROMIC GUT
101-2787	C-6	27"	4-0	CHROMIC GUT
101-8824	C-3	18"	5-0	CHROMIC GUT
900-4417	C-3	18"	5-0	PGA - UNDYED BRAIDED
900-4418	C-3	18"	4-0	PGA - UNDYED BRAIDED
900-4421	PC-31	18"	3-0	PGA - UNDYED BRAIDED
900-4419	PC-31	18"	4-0	PGA - UNDYED BRAIDED
900-4420	C-6	18"	4-0	PGA - UNDYED BRAIDED
900-7462	C-6	27"	3-0	PGA - UNDYED BRAIDED
900-7463	C-6	27"	4-0	PGA - UNDYED BRAIDED
900-7468	C-3	18"	6-0	PGA - UNDYED BRAIDED
900-7477	PC-30	18"	4-0	PGA - UNDYED BRAIDED
900-7465	C-3	18"	4-0	PGA-PCL-MONO
900-7467	C-3	18"	5-0	PGA-PCL-MONO
900-7474	PC-31	18"	3-0	PGA-PCL-MONO
900-7488	DSM-19	18"	4-0	PGA-PCL-MONO

**SUTURES
NON-ABSORBABLE**

All Sutures..... 12/Box

ITEM CODE	NEEDLE	LENGTH	SIZE	MATERIAL
101-6053	C-22	18"	6-0	BLACK NYLON
101-7137	C-6	18"	5-0	BLACK NYLON
101-1209	C-6	18"	4-0	BLACK NYLON
101-3728	C-7	18"	3-0	BLACK NYLON
900-7459	C-17	10"	5-0	BLACK NYLON
900-7460	C-17	10"	6-0	BLACK NYLON
900-7461	C-17	10"	4-0	BLACK NYLON
900-7464	C-6	10"	4-0	BLACK NYLON
900-7466	C-3	18"	4-0	BLACK NYLON
900-7470	C-3	18"	5-0	BLACK NYLON
900-7472	PC-31	18"	5-0	BLACK NYLON
900-7473	C-3	18"	6-0	BLACK NYLON
900-7476	C-8	18"	2-0	BLACK NYLON
900-7482	C-1	18"	6-0	BLACK NYLON
900-7487	PC-31	18"	4-0	BLACK NYLON
101-4378	C-22	18"	6-0	BLUE POLYPROPYLENE
101-6409	C-6	18"	5-0	BLUE POLYPROPYLENE
101-7477	C-6	18"	4-0	BLUE POLYPROPYLENE
101-7700	C-7	18"	3-0	BLUE POLYPROPYLENE
900-7478	PC-31	18"	4-0	BLUE POLYPROPYLENE
900-7479	C-3	18"	6-0	BLUE POLYPROPYLENE
900-7481	DS-18	18"	3-0	BLUE POLYPROPYLENE
900-7483	C-1	18"	6-0	BLUE POLYPROPYLENE
900-7484	C-3	18"	5-0	BLUE POLYPROPYLENE
900-7486	C-3	18"	4-0	BLUE POLYPROPYLENE
100-7762	C-43	18"	3-0	BLACK BRAIDED SILK
100-6320	C-9	18"	3-0	BLACK BRAIDED SILK
100-7772	C-6	18"	3-0	BLACK BRAIDED SILK
100-6830	C-6	18"	4-0	BLACK BRAIDED SILK
101-3207	C-6	18"	5-0	BLACK BRAIDED SILK
101-2590	C-3	18"	4-0	BLACK BRAIDED SILK
101-4148	C-3	18"	5-0	BLACK BRAIDED SILK
101-2049	C-7	18"	3-0	BLACK BRAIDED SILK
101-2636	C-22	18"	6-0	BLACK BRAIDED SILK
100-7773	D-14	18"	5-0	GREEN BRAIDED POLYESTER

Podiatry Innovative products you can trust!




"PBS" Diamond Burs
Nail Reduction and
Tissue Debridement

Professional treatment system for hyperkeratotic tissue and hypertrophic nails.



Barrel • 378-8345
(To reduce thick, deformed nails.)



Pear • 378-7071 ...
(To smooth the rough edges of nails.)



Bud • 378-8342 ...
(The most popular pattern for nail reduction.)

179-7955
Diamond Bur System Kit ...

A convenient kit of four nail and callus skin burs.
(Includes one each: Bud, Short, Taper, Umbrella.)




Short • 325-4739
(For better control in treating smaller digits.)



Olive • 101-3186
(Similar to the PEAR and BUD.)



Taper • 378-9386
(For the treatment of incurvated nails or the nail groove.)



Umbrella • 378-5542 ...
(For the reduction of hyperkeratoses, tylomas and heel callus.)



Small Olive • 104-9402 .
(For working smaller areas.)



Small Cylinder • 893-0006
(For gross reduction of large nail areas.)

Lumicain®
Hemostatic Agent

378-4603
Lumicain®



Salinocaine™
For treatment of verruca

496-3421
Salinocaine



Accu-Mold®
Moldable Silicone Compound



For custom digital appliances
378-3225
Accu-Mold
200 Gram Kit (Yields approximately 15 devices.)



SWISS CILIA AND SUTURE/SPLINTER FORCEPS 3¾"

German, Precision Fitted, Diagonal Jaws, Stainless Steel
Great for hair and suture removal.
(102-3436)ea



LISTER BANDAGE SCISSORS 4½"
German, Extra Fine
(104-8240)ea



KNOWLES BANDAGE SCISSORS 5½"
German
(104-5504)ea



LISTER BANDAGE SCISSORS 7¼"
German, Extra Fine
(101-5267)ea

HENRY SCHEIN®
BRAND PRODUCTS

LISTER BANDAGE SCISSORS 7½"
(104-7364)ea



KELLY FORCEPS 5½"
Curved
(100-7878)ea



OLSEN-HEGAR NEEDLE HOLDER 5½"
German, Suture/Scissors Combo, Serrated Jaws
(100-1395)ea



IRIS SCISSORS 3½"
German, Straight, Delicate Pattern
(102-5248)ea



TAPE 'N TELL

Eight different-colored self-adhering strips for instrument identification. Autoclavable.

Kit
(100-8919)ea

Contains: 1 of ea tape: orange, blue, brown, green, white, black, red & yellow.
9½" Refill Rollsea

Specify:
Orange(294-1208)
Blue(294-3969)
Brown(294-4335)
Green(294-4712)
White(294-4979)
Black(294-5260)
Red(294-5400)
Yellow(294-5735)



SELF-RETAINING RETRACTOR 4"
Sharp, Blunt Points
(104-5534)ea

Disposable Needles and Syringes



HENRY SCHEIN® SYRINGES WITHOUT NEEDLES

1 cc (900-4480)	100/box
3 cc (900-4475)	100/box
5 cc (900-4477)	100/box
10 cc (900-4476)	100/box
20 cc (112-6151)	50/box
30 cc (112-6153)	50/box
60 cc (900-4479)	50/box

HENRY SCHEIN®
BRAND PRODUCTS

HENRY SCHEIN® DISPOSABLE NEEDLES

30 ga x 1" (112-6149)	100/box
30 ga x 1/2" (900-4467)	100/box
25 ga x 5/8" (900-4468)	100/box
25 ga x 1" (900-4470)	100/box
25 ga x 1 1/2" (900-4472)	100/box
23 ga x 1" (900-4473)	100/box
18 ga x 1" (900-4469)	100/box
18 ga x 1 1/2" (900-4471)	100/box

HENRY SCHEIN® LUER LOCK 3-CC DISPOSABLE SYRINGES WITH NEEDLES

Syringes	100/box
Specify:	
25 ga x 5/8"	(900-4458)
25 ga x 1"	(900-4461)
23 ga x 1"	(900-4460)
21 ga x 1"	(900-4462)

Sharps Collectors



- Ergonomic handle for the 1.0-quart/1.5-quart/5-quart and 2-gallon and 3 gallon collectors
- Anti-defluxion wing sizes prevent outflow of waste
- Maximum filling level indicator
- No emission of any harmful gases during the incineration process
- Simple and practical snap-type closure
- Easy maneuvering for temporary opening with just one hand
- Wide opening for the disposal of large amounts of sharps waste
- Mechanism for the removal of needles on some units
- 4-quart/8-quart/18-gallon models are ideal for the disposal of long needles

Red	ea
(900-4980)	1 qt
(900-4981)	1.5 qt
(900-4982)	4 qt
(900-4983)	5 qt
(900-4984)2 gal
(900-4985)3 gal
(900-4986)8 gal
(900-7032)	18 gal (7 pk)

HENRY SCHEIN®
BRAND PRODUCTS



**CRITERION®
PC POWDER-FREE
LATEX EXAM GLOVES**

Non-sterile, Polymer Coated, White
Textured • Nonchlorinated
Sets the standard in quality, powder-free gloves. Nonchlorinated, synthetic polymer coating eases donning. Textured finish for enhanced grip. Excellent tactile sensitivity.
Gloves 100/box

- Specify:
- X-small(102-5418)
 - Small(102-5419)
 - Medium.....(102-5421)
 - Large(102-6730)
 - X-large (90/box).....(102-5422)



**CRITERION®
POWDER-FREE
NITRILE EXAM
GLOVES**

Non-sterile, Blue
Fully textured nitrile glove that offers excellent tactile sensitivity. Soft and stretchy, this nonchlorinated glove is the ideal alternative to latex.
Fully Textured • Nonchlorinated
Gloves 100/box

- Specify:
- X-small(111-8533)
 - Small(111-8535)
 - Medium.....(111-8536)
 - Large(111-8537)
 - X-large (90/box).....(111-8539)



**CRITERION®
POWDER-FREE
VINYL EXAM
GLOVES**

Non-sterile
Contains no natural latex, providing an excellent solution for those at risk of latex sensitivities.

- Supple touch for maximum tactile sensitivity
- Powder-free to eliminate glove-powder aerosolization

Smooth
Gloves 100/box

- Specify:
- X-small(102-5721)
 - Small(102-6721)
 - Medium.....(102-5338)
 - Large(102-7319)
 - X-large (90/box).....(102-5377)

HENRY SCHEIN®

BRAND PRODUCTS

**CRITERION®
POWDER-FREE
LATEX
SURGEONS'
GLOVES**



Sterile
Henry Schein's Criterion brand of latex surgical gloves helps eliminate powder-related complications and irritation, yet still provides outstanding quality and performance you would expect from a premium surgical glove.
• Supersoft latex film and anatomical design provide unprecedented comfort to reduce fatigue, even during the longest procedures

Textured • Chlorinated • Certified Protein
Claim: None

Gloves.....50 pair/box

- Specify:
- 5½(104-6125)
 - 6(104-5073)
 - 6½(104-9218)
 - 7(104-3502)
 - 7½(104-3503)
 - 8(104-3563)
 - 8½(104-7315)
 - 9(104-7227)

HAND SANITIZER

- Alcohol-based waterless hand gel
- Kills 99.9% of most common disease-causing germs spread by the hands
- Rinse-free; no water required
- Fresh clean scent and soothing emollients leave hands soft and clean



8-oz Pump Bottle
(900-4443).....ea

16-oz Pump Bottle
(900-4442).....ea

**FOAMING
HAND
SANITIZER**

Alcohol-free, waterless foaming hand sanitizer. Special formulation contains emollients and a fresh pleasant Mountain Meadow fragrance, leaving hands feeling soft and clean after use, without any stickiness or drying out of the skin. Rinse-free; no water required. Kills 99.9% of the most common disease-causing germs spread by hands. Contains benzalkonium chloride (active ingredient) to kill germs.



18.6-oz Pump Bottle
(900-4444).....ea

**ANTIBACTERIAL
LIQUID HAND SOAP**

Luxurious Liquid Hand Soap for Sinks and Dispensers
Provides a luxurious feel, leaving hands smooth and soft. Mild, gentle, and moisturizing. Pleasant Cucumber and Aloe fragrance. Active ingredient is 0.2% triclosan.



- 8-oz Pump Bottle
(900-4439).....ea
- 16-oz Pump Bottle
(900-4438).....ea
- 1-gal Refill Bottle
(900-4440).....ea

**MAXICLENS®
ANTIMICROBIAL/
ANTISEPTIC
SKIN CLEANSER**

Formulated with detergents, moisturizers, skin protectors, and humectants that cleanse without drying. Fast-acting formula that is proven effective at 15 seconds against a broad spectrum of bacteria. Effective for up to 6 hours after use. Active ingredient: 4% chlorhexidine gluconate.



- 1 Quart
(431-0032).....ea
- 1 Gallon
(431-0099).....ea
- Hand Pump for Quart Bottle
(431-1478).....ea

QUICK WIPES

Nonallergenic, premoistened, ready-to-use hand towelettes. Pleasant fragrance with lanolin, alcohol, and aloe. Refreshing aid for the office, home, or car. Use for easy, gentle cleansing of hands and face.



Pop-Up Container, 9" x 5½"
(104-4176).....150/cont

HAND LOTION

- Luxurious hand cream
- Soothing to skin
- Contains emollients to leave hands feeling soft
- Pleasant Cucumber and Aloe fragrance



8-oz Pump Bottle
(900-4446).....ea

Gebauer's Ethyl Chloride®


Topical Anesthetic Skin Refrigerant

A Proven Podiatric Solution

Is your current generic cold spray safe and proven like Gebauer's Ethyl Chloride?

Podiatrists have trusted Gebauer's Ethyl Chloride for more than 100 years. It is FDA cleared and clinically supported, so you can be confident you are offering your patients a safe and effective topical anesthetic. Gebauer's Ethyl Chloride instantly and temporarily controls the pain associated with needle procedures and minor surgical procedures, like:

- Podiatric Blocks
- Lidocaine Injections
- Ingrown Toenails
- Incision and Drainage of small abscesses
- Lancing Boils
- Skin Tag Removal



Important Risk and Safety Information for Gebauer's Ethyl Chloride:

Published clinical trial results support the use in children 3 years of age and older. Ethyl Chloride is FLAMMABLE and should never be used in the presence of an open flame or electrical cautery equipment. Inhalation should be avoided as it may produce narcotic and general anesthetic effects, and may produce deep anesthesia, fatal coma or cardiac arrest. Do not spray in eyes. Over spraying may cause frostbite. Freezing may alter skin pigmentation. The thawing process may be painful, and freezing may lower resistance to infection and delay healing. Cutaneous sensitization may occur, but appears to be extremely rare. Long-term exposure may cause liver or kidney damage. Rx Only.



View a lecture and get CME credit

Via PRESENTS e-Learning Systems (registration required)

"Increasing Podiatric Patient Satisfaction Through the Use of Topical Anesthetics"
by Dr. Scott Hirsh, DPM, DABPS

Learn different techniques to increase patient satisfaction by using topical anesthetics to reduce pain and also to reduce patient stress levels.

www.Gebauer.com/PatientSatisfaction



Gebauer Company

The Most Trusted Name In Skin Refrigerants
For Over 100 years!®

www.Gebauer.com

Products Made in the U.S.A.

4444 East 153rd Street • Cleveland, Ohio 44128
1-800-321-9348

Betamethasone Sodium Phosphate and Betamethasone Acetate

Compare to Celestone® Soluspan®
6 mg/mL, 5-mL MDV
(114-8668)ea



Celestone® Soluspan® Injectable Suspension

6 mg/mL, 5-mL MDV
(516-0266)ea



Depo-Medrol® Injectables With Preservative

20 mg/mL, 5-mL MDV
(908-3069)ea



Depo-Medrol® Injectables Methylprednisolone Acetate

40 mg/mL, 10-mL MDV
(908-1188)ea
80 mg/mL, 5-mL MDV
(908-3787)ea



Dexamethasone Sodium Phosphate Injectables

Compare to Decadron®
4 mg/mL, 30-mL MDV
(102-4486)ea



Gebauer's Ethyl Chloride®

Topical Anesthetic Skin Refrigerant,
3½ oz Glass Bottle with Nozzle
Fine-Stream Spray
(761-1189)4/box
Medium-Stream Spray
(761-9264)4/box



Formadon

10% Formaldehyde USP Strength
Drying agent and antiperspirant used in the treatment of excessive conditions of odor and wetness. Also treats mosaic warts.

2 oz
(101-8108)ea

Gentamicin Cream

0.1%, 15-g Tube
(116-7311)ea

Gentamicin Cream

Compare to Garamycin®
0.1%, 30-g Tube
(108-7056)ea



Gordal Creme 20% Urea

2½ oz
(812-9583)ea
4 oz
(812-9924)ea

Kenalog®-10 Injectable

Triamcinolone Acetonide
10 mg/mL, 5-mL MDV
(196-0545)ea



Kenalog®-40 Injectables

Triamcinolone Acetonide
40 mg/mL, 5-mL MDV
(196-1737)ea
40 mg/mL, 10-mL MDV
(196-9429)ea



Lidocaine HCl Injectable

Compare to Xylocaine®
50-mL MDVs
1%
(104-6817)25/box
2%
(104-9843)25/box



PREMIER

Lumicain® Solution

For rapid control of minor hemorrhages. Because it is colorless, Lumicain provides a clear vision of the operating area.

60 cc
(378-4603)ea



Marcaïne Injectables

0.25%, 10-mL SDV, Preservative-free
(631-2127)10/box
0.5%, 50-mL MDV
(631-2615)ea



Marcaïne® with Epinephrine Injectables

0.5%, 50-mL MDV
(258-9815)ea



Mycomist

Fungicidal, sanitizing deodorant spray for shoes and boots. Utilizes the triple action of chlorophyll, formalin, and benzalkonium chloride.

4-oz Bottle
(108-2471)ea



Silvadene® Cream

20-g Tube
(719-8272)ea



Silver Nitrate Applicator

- 6" plastic applicator stick
 - 75% silver nitrate, 25% potassium nitrate
 - Automated production provides consistent silver nitrate on applicator
 - Packaged in a secure vial to ensure moisture control
 - Prescription item
- (277-0017)100/vial



Tri-Chlor

Trichloroacetic acid in a clear liquid solution. It is used as a cauterant to treat condyloma.

80%, 15-cc Bottle
(812-5871)ea

NEW

Introducing the
Every Which Way.
Biofreeze® 360° Spray

Best. Biofreeze. Ever.
Now in gel, roll-on, and spray
BIOFREEZE
COLD THERAPY PAIN RELIEF



UNIFORM SPRAY PATTERN > for optimal coverage
CONTINUOUS SPRAY > much easier than repeated pumping
360° SPRAY > best spray technology... works from any angle

Biofreeze® 360° Spray
4 oz
(601-0871)ea **\$7.69**

biofreeze.com | 800.246.3733 | Like ice, only better

Biofreeze is a registered trademark of The Hygenic Corporation. Unauthorized use is strictly prohibited.
©2012 The Hygenic Corporation. All rights reserved.



ExSept Plus Skin and Wound Cleanser

Electrolytically Produced Sodium Hypochlorite 0.114% Sodium Hypochlorite Intended for the cleansing of skin debris and foreign material from wounds.



100-mL Bottle (602-9874)24/case
200-mL Spray (602-7365)ea
250-mL Bottle (602-1426)ea
500-mL Bottle (602-1218)12/case



Foot Miracle Cream

This unique, deep-penetrating formula softens calloused rough areas, bringing comfort and vitality to heels and feet. Also recommended for elbows, knees and any other rough-skin area. Never sticky or greasy. Leaves skin feeling soft.

6 oz (614-9342)ea

Urea Care

6 oz (180-2216)ea



Gordon's #5 Foot Powder Spray

4 oz (904-3283)ea



Emollia Lotion

4 oz (763-7932)ea
8 oz (366-7897)ea



Gormel Creme 20% Urea

2½ oz (812-9583)ea
4 oz (812-9924)ea



Aloe Grande Creme

Fine moisturizing skin cream promotes postoperative healing. Contains vitamin A 100,000 IU per ounce, vitamin E 1500 IU per ounce and aloe vera gel.

4 oz (987-2536)ea



DaniPro Antifungal Nail Polish

- Infused with an FDA-approved antifungal
- Organic and 100% natural
- Safe and effective
- 4Free formulation:
 - No formaldehyde
 - No DBP
 - No toluene
 - No animal testing

Antifungal Nail Polish
½-oz Bottleea

Specify:
Brown ("Someone to Love")(115-5997)
Clear Base ("Be Free")(115-5996)
Clear Top ("Peace")(115-5994)
Deep Red ("Always Be True")(115-6000)
Mocha ("Babe It's You")(115-5999)
Pink ("Love Is All")(115-6002)
Red ("First Kiss")(115-6001)
Wine ("Say You Love Me")(115-6007)
Lilac ("Stop and Smell the Flowers")(116-6657)
Nude ("Nothing to Hide")(116-6644)
Steel ("Show Your Strength")(116-6667)
Pure Pink ("My Girl") NEW!(125-0133)
True Blue NEW!(365-0293)
Cocoa Cabernet NEW!(365-0292)

Free Display! Call your Rep for Details!

Offer ends March 31, 2013

MENTION PROMO CODE 4F
to ensure promotional pricing.

1-800-323-5110

8am-9pm (et)

Fax:

1-800-524-4969

24 Hours

www.henryschein.com/footandankle

Due to manufacturer price increases, prices are subject to change without notification. Offer valid on our full array of products, excluding flu vaccine and Synagis®. Upon any request by Medicare, Medicaid, Tricare or other payer, you must report the value of any benefit received from a discount program (e.g., points redemptions, gifts or other special awards). ©2013 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. No adjustments from prior sales. Offer not to be combined with other promotions and/or special contract pricing. Void where prohibited. Local regulatory requirements may apply to use or installation of certain products. It is your responsibility to understand and comply with any such requirements prior to purchase, use or installation of products. †For informational purposes only. Source: ReimbursementCodes.com and/or www.cms.gov. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.



Printed on recycled paper.
Please recycle it again.

13MS7514

HENRY SCHEIN®

Henry Schein, Inc. • 135 Duryea Road • Melville, NY 11747

PRSR STD
U.S. POSTAGE

PAID

HENRY SCHEIN, INC.



Which Rewards Program Will You Choose?

Cashback Rebate Rewards

- Earn **2% CASH BACK** on qualifying Henry Schein transactions up to \$100,000 a year†
- Earn **UP TO 2% CASH BACK** on other qualifying transactions†
- Get **AUTOMATIC REDEMPTIONS** in \$50 increments, posted to your account†



Maximum Rewards®

- Earn **1.5 POINTS** per \$1 in qualifying Henry Schein transactions*
- Earn **1 POINT** per \$1 in other qualifying transactions*
- Redeem for **YOUR CHOICE** of exciting travel, merchandise, or gift cards

Both Options Offer:

- **Extended Protection Program for Equipment Purchases**
- **Full Travel Benefits • No Annual Fee**

HENRY SCHEIN®

FINANCIAL SERVICES

Business Solutions for Healthcare Practitioners

1-866-398-9296 or
www.henryschein.com/creditcard

The issuer of the card reserves the right to change the credit card account terms should your credit profile change or for other reasons described in the Cardmember Agreement. Cards are issued by First Bankcard, a division of First National Bank of Omaha.

†You will receive 2% cash back on up to \$100,000 annually on qualifying Henry Schein transactions and 0.50% on qualifying Henry Schein transactions above \$100,000.

*Complete details, including restrictions, limitations and exclusions will be available when you become a Cardmember.

Henry Schein Financial Services is not a bank, does not represent itself as such, and does not conduct banking activities.
© 2013 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors.