

BLU-BITE® HP
FAST & SUPER-FAST SET

**EXPERIENCE THE
ACCURACY.**

NEW
**OUR MOST
RIGID
FORMULA**

**CAPTURE CONSISTENT
AND RELIABLE
BITE IMPRESSIONS**

with Blu Bite® HP registration material in fast or super-fast setting times. Used for closed, open, and double-bite registrations.



For more information, see product Instructions for Use.

For our HENRY SCHEIN BRAND PRODUCTS CATALOG, as well as additional product information, visit WWW.HENRYSCHEINDENTAL.COM/HSBRAND.



VALUE & QUALITY



RELIABLE & TRUSTWORTHY



PRODUCT SELECTION

BLU-BITE® HP

FAST & SUPER-FAST SET



NEW

OUR MOST RIGID FORMULA

Advantages

Provides low resistance to closure and minimizes deviation of the mandible for accuracy and reliability

Increases patient's comfort with pleasant scents: Chocolate, Vanilla, and Mint

Fast-setting, easy-to-extrude and trim



Henry Schein® Blu-Bite® HP

COMPLETE KIT 2/Box + 3 T-Mix Tips

Scent	Rigidity	Set Time	Item Code	Color
Mint	Shore D > 40	0:55	570-2967	Teal
Vanilla	Shore A 90	0:55	570-2966	White
Chocolate	Shore A 90	1:30	570-2965	Brown

COMPATIBLE MIXING TIP (48/Bg) T-Mix 112-6862

For more information, see product Instructions for Use.

For our HENRY SCHEIN BRAND PRODUCTS CATALOG, as well as additional product information, visit WWW.HENRYSCHEINDENTAL.COM/HSBRAND.



VALUE & QUALITY



RELIABLE & TRUSTWORTHY



PRODUCT SELECTION