

CENTRIX
Absolute



PRODUCT REVIEW SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Ease of dispensing/placement
with intraoral tips **4.6**

Viscosity **4.2**

Working time **4.4**

Curing/setting time **4.2**

Esthetics—shade match &
shade choices available **4.0**

Ease of cleanup **4.0**

Lack of postop sensitivity **4.6**

SECTION A
AVERAGE **4.3**

SECTION B
OVERALL SATISFACTION **4.1**

takeaways

- Available in 5-g syringes with automatic mixing tips and snap-on intraoral tips
- Working time of 1 minute 30 seconds and a setting time of 4 minutes
- 3 shades (opaque, translucent, and universal) for a natural-looking restoration
- High flexural and compressive strengths ensure superior adhesion to dentin and ceramic
- Easy cleanup saves chair time

FINAL SCORE
AVERAGE OF SECTIONS A AND B

4.2