

DENTSPLY SIRONA  
**SmartLite Pro**



**PRODUCT REVIEW** SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Ergonomics **4.5**

Curing time/speed of cure **4.6**

Access to and visibility of oral cavity **4.0**

Complete coverage of area being cured **4.6**

Portability **5.0**

Battery life **4.6**

Design of charging base **4.2**

SECTION A  
AVERAGE **4.5**

SECTION B  
OVERALL SATISFACTION **4.4**

## takeaways

- State-of-the-art LED optics provide quality cures for reliable outcomes
- Lightweight design with a compact pen-style grip
- Output of 1,250 mW/cm<sup>2</sup>, curing width of 10 mm, curing depth of 8 mm
- Interchangeable, 360° rotational tips for a variety of clinical indications
- Multifunctional charging base with built-in radiometer, extra battery, and storage for extra tips

**FINAL SCORE**  
AVERAGE OF SECTIONS A AND B

**4.5**