

DIRECTA DENTAL

Luxator P-series



PRODUCT REVIEW SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Cutting ability **4.7**

Ability to extract without trauma **4.5**

Ease of access to oral cavity (due to small size) **4.6**

Ergonomics—grip, weight, pen-like design **4.6**

Sharpness of blade **4.8**

Reduced friction due to titanium tip **4.6**

Selection of sizes **4.4**

SECTION A
AVERAGE **4.6**

SECTION B
OVERALL SATISFACTION **4.7**

takeaways

- Thin, sharp periotome tips efficiently cut the periodontal ligament before extraction
- Lightweight, pen-like grip for improved ergonomics and tactility
- Blades coated with titanium nitride for durability; no sharpening needed
- Four blade configurations for any clinical situation—straight, inverted curve, contra angle, and dual edged
- Instrument design enables fast, atraumatic extractions

FINAL SCORE
AVERAGE OF SECTIONS A AND B

4.7