

# HENRY SCHEIN DENTAL ADAPTA-Post



## PRODUCT REVIEW SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Ease of post insertion/  
placement **4.5**

Venting of cement **4.5**

Fit between post and coronal section **4.4**

Sharpness of threads—allows for smooth cuts and does not compress dentin **4.4**

Radiopacity of post **4.6**

Retention of core buildup material to post (initial retention in office day of procedure) **4.6**

SECTION A  
AVERAGE **4.5**

SECTION B  
OVERALL SATISFACTION **4.4**

## takeaways

- Tapered tip creates smooth transition between post and root canal
- Split-shank design transmits insertion stress toward post instead of root
- Sharp threads on post cut smoothly into dentin
- Available in titanium and stainless steel in sizes 0, 1, 2, and 3
- Wide coronal shank ensures a proper fit between post and coronal section

**FINAL SCORE**  
AVERAGE OF SECTIONS A AND B

**4.5**