

MICROCOPY

# NeoDiamond Z-Class Adjusters



**PRODUCT REVIEW** SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Level of chatter **4.4**

Ability to make accurate  
& controlled adjustments **4.4**

Polishing time required **4.2**

Final finished surface after  
using adjuster and polisher **4.2**

Packaging & labeling **4.7**

Range/variety of shapes **4.5**

SECTION A  
AVERAGE **4.4**

SECTION B  
OVERALL SATISFACTION **4.4**

## takeaways

- Z-Class diamonds are designed specifically for cutting and adjusting zirconia restorations
- Sterile, single-patient-use, zirconia-grade diamonds
- Produces a smooth cut at a consistent depth without scarring
- Because of the smooth cutting during adjustment, polishing time is reduced
- Adjusters come in 7 shapes—round end taper, ball, teardrop, needle, egg, flame, pointed cone

**FINAL SCORE**  
AVERAGE OF SECTIONS A AND B

**4.4**