

TOKUYAMA DENTAL

# Estelite Bulk Fill Flow



## PRODUCT REVIEW SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Delivery system **4.6**

Procedure time **4.6**

Cavity adaptation **4.6**

Overall ease of use and simplicity of procedure **4.4**

Polishing and glossiness **4.3**

Shade matching and final esthetics **3.6**

SECTION A  
AVERAGE **4.4**

SECTION B  
OVERALL SATISFACTION **4.1**

## takeaways

- Quick procedure with no capping layer and no incremental layering necessary
- Good cavity adaptation, flows into cracks, excellent marginal seal
- 5 shades with superior color matching ability
- Quickly cures to depths of 4 mm
- Low shrinkage stress, high compressive strength

**FINAL SCORE**  
AVERAGE OF SECTIONS A AND B

**4.3**