


VP MIX™ HP MAX HYDRO
IMPRESSION MATERIAL

**ACHIEVE
MAXIMUM
DETAIL WITH
INSTANT
HYDROPHILICITY.**

Leverage our highest level of hydrophilicity to make a great impression the first time and take the first step in a high-quality, low-discomfort restoration.



For more information, see product Instructions for Use.

Order today at www.henryscheindental.com/HSBImpression



VALUE & QUALITY



RELIABLE & TRUSTWORTHY




PRODUCT SELECTION

VP MIX™ HP MAX HYDRO IMPRESSION MATERIAL

CONTACT TIME (Seconds)

0
1
3
5
10
20
30



Advantages

A low-contact angle delivers maximum accuracy of a patient's dentition

High-contrast colors and saturation

Maximum tear strength and elasticity maintains the impression's integrity during removal

Pleasant scents reduce patient gag reflex and discomfort

As shown on the left, water droplets on VP Mix™ HP Max Hydro dropped to 0° within 30 seconds. By definition, a contact angle of 0° is considered super hydrophilic.

Unflavored



HP MAX HYDRO, INTRO KIT

2 Cartridges (48ml each) + 3 T-Mix Tips

SET TYPE & TIME	VISCOSITY	COLOR	ITEM NUMBER
Fast Set (2:15)	Extra Light Body	Orange	(570-2988)
	Light Body	Green	(570-2990)
	Regular Body	Lilac	(570-2992)
	Heavy Body	Purple	(570-2994)
	Monophase	Blue	(570-2996)
Regular Set (4:30)	Extra Light Body	Orange	(570-2989)
	Light Body	Green	(570-2991)
	Regular Body	Lilac	(570-2993)
	Heavy Body	Purple	(570-2995)
	Monophase	Blue	(570-2997)

360 MAX HYDRO, JUMBO

(360ml) + 5 Dynamic Mixing Tips + 1 Bayonet Ring

Regular Set	Heavy Body	Purple	(570-3007)
	Monophase	Blue	(570-3006)

HP MAX HYDRO, BULK REFILL

8 Cartridges (48ml each, no tips included)

SET TYPE & TIME	VISCOSITY	COLOR	ITEM NUMBER
Fast Set (2:15)	Extra Light Body	Orange	(570-3033)
	Light Body	Green	(570-3035)
	Regular Body	Lilac	(570-3037)
	Heavy Body	Purple	(570-3039)
	Monophase	Blue	(570-3041)
Regular Set (4:30)	Extra Light Body	Orange	(570-3034)
	Light Body	Green	(570-3036)
	Regular Body	Lilac	(570-3038)
	Heavy Body	Purple	(570-3040)
	Monophase	Blue	(570-3042)

PUTTY MAX HYDRO, 2-PACK

(250ml each) + 2 Measuring Spoons

Regular Set	Heavy Body	Purple	(570-3001)
Fast Set	Heavy Body	Blue	(570-3000)

For more information, see product Instructions for Use.

Order today at www.henryscheindental.com/HSBImpression



VALUE & QUALITY



RELIABLE & TRUSTWORTHY



PRODUCT SELECTION

VP MIX™ HP MAX HYDRO IMPRESSION MATERIAL

Berry



HP MAX HYDRO, INTRO KIT

2 Cartridges (48ml each) + 3 T-Mix Tips

SET TYPE & TIME	VISCOSITY	COLOR	ITEM NUMBER
Fast Set (2:15)	Extra Light Body	Orange	(570-2968)
	Light Body	Green	(570-2970)
	Regular Body	Lilac	(570-2976)
	Heavy Body	Purple	(570-2972)
	Monophase	Blue	(570-2974)
Regular Set (4:30)	Extra Light Body	Orange	(570-2969)
	Light Body	Green	(570-2971)
	Regular Body	Lilac	(570-2977)
	Heavy Body	Purple	(570-2973)
	Monophase	Blue	(570-2975)

360 MAX HYDRO, JUMBO

(360ml) + 5 Dynamic Mixing Tips + 1 Bayonet Ring

Regular Set	Heavy Body	Purple	(570-3004)
	Monophase	Blue	(570-3005)

HP MAX HYDRO, BULK REFILL

8 Cartridges (48ml each, no tips included)

SET TYPE & TIME	VISCOSITY	COLOR	ITEM NUMBER
Fast Set (2:15)	Extra Light Body	Orange	(570-3023)
	Light Body	Green	(570-3025)
	Regular Body	Lilac	(570-3029)
	Heavy Body	Purple	(570-3027)
	Monophase	Blue	(570-3031)
Regular Set (4:30)	Extra Light Body	Orange	(570-3024)
	Light Body	Green	(570-3026)
	Regular Body	Lilac	(570-3030)
	Heavy Body	Purple	(570-3028)
	Monophase	Blue	(570-3032)

PUTTY MAX HYDRO, 2-PACK

(250ml each) + 2 Measuring Spoons

Regular Set	Heavy Body	Purple	(570-2998)
Fast Set	Heavy Body	Blue	(570-2999)

Mint



HP MAX HYDRO, INTRO KIT

2 Cartridges (48ml each) + 3 T-Mix Tips

SET TYPE & TIME	VISCOSITY	COLOR	ITEM NUMBER
Fast Set (2:15)	Extra Light Body	Orange	(570-3008)
	Light Body	Green	(570-3010)
	Regular Body	Lilac	(570-3014)
	Heavy Body	Purple	(570-3012)
	Monophase	Blue	(570-3016)
Regular Set (4:30)	Extra Light Body	Orange	(570-3009)
	Light Body	Green	(570-3011)
	Regular Body	Lilac	(570-3015)
	Heavy Body	Purple	(570-3013)
	Monophase	Blue	(570-3017)

PUTTY MAX HYDRO, 2-PACK

(250ml each) + 2 Measuring Spoons

Regular Set	Heavy Body	Purple	(570-3002)
Fast Set	Heavy Body	Blue	(570-3003)

HP MAX HYDRO, BULK REFILL

8 Cartridges (48ml each, no tips included)

SET TYPE & TIME	VISCOSITY	COLOR	ITEM NUMBER
Fast Set (2:15)	Extra Light Body	Orange	(570-3043)
	Light Body	Green	(570-3045)
	Regular Body	Lilac	(570-3047)
	Heavy Body	Purple	(570-3049)
	Monophase	Blue	(570-3051)
Regular Set (4:30)	Extra Light Body	Orange	(570-3044)
	Light Body	Green	(570-3046)
	Regular Body	Lilac	(570-3048)
	Heavy Body	Purple	(570-3050)
	Monophase	Blue	(570-3052)

For more information, see product Instructions for Use.

Order today at www.henryscheindental.com/HSBImpression



VALUE & QUALITY



RELIABLE & TRUSTWORTHY



PRODUCT SELECTION

VP MIX™ HP MAX HYDRO



T-Mix Tips



Helical Mixing Tips



Intraoral Tips



Dynamic Mixing Tips



VP Mix Bayonet Ring

Mixing Tips Reordering Guide

FOR USE WITH VP MIX HP MAX HYDRO

VISCOSITIES	DESCRIPTION	TIP COLOR / SIZE	QUANTITY	ITEM NUMBER
Extra Light, Light & Regular Body	T-Mix Tips	Yellow/4.2mm	48/bag	(112-6863)
	Helical Mixing Tips	Yellow/4.2mm	48/bag	(900-0764)
	Intraoral Mixing Tips	Yellow	100/bag	(900-4587)
Monophase & Heavy Body	T-Mix Tips	Green/6.5 mm	48/bag	(112-6862)
	Helical Mixing Tips	Green/6.5 mm	48/bag	(900-0766)
	Intraoral Mixing Tips	Clear	100/bag	(112-6704)

FOR USE WITH VP MIX 360 MAX HYDRO, JUMBO

Dynamic Mixing Tips*	Yellow	50/bag	(900-4598)
VP Mix Bayonet Ring	Yellow	2/pk	(900-4599)

*380ml mixing tips are compatible with 360 MAX HYDRO, JUMBO cartridges.

For more information, see product Instructions for Use.

Order today at www.henryscheindental.com/HSBImpression



VALUE & QUALITY



RELIABLE & TRUSTWORTHY



PRODUCT SELECTION