

ZIRC

Performance Tool Organizer



PRODUCT REVIEW SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Overall fit of instruments in the organizer **3.8**

Fit of organizer in ultrasonic, autoclave, or Statim **4.0**

Protection of dental tools and dental professional **4.3**

Color options available for color coding by procedure **4.0**

Overall design of organizer **4.3**

Water flow **4.3**

Overall efficiency **3.7**

SECTION A AVERAGE **4.1**

SECTION B OVERALL SATISFACTION **3.8**

takeaways

- Ideal for large and specialty instruments
- Customizable, removable feet to hold instruments in place
- Compact, lightweight design to fit most sterilizers
- Attached cover rolls back easily
- Many colors available and dimensions are 8-3/16" x 4-1/4" x 1-7/32"

FINAL SCORE
AVERAGE OF SECTIONS A AND B

4.0

Dental **PRODUCT SHOPPER**

RECOMMENDED
PRODUCT