



# POINT-OF-CARE



## OSOM® Ultra Plus Flu A&B Test

### STRONGER CLINICAL PERFORMANCE TAKES LATERAL FLOW TESTING TO THE NEXT LEVEL

#### Features and Benefits


- **Accurate** – High Sensitivity / Specificity
- **CLIA-waived** – No need for instrument
- **Fast and simple workflow** – Pre-filled Vials, Minimal hands on time, Results in 10 minutes
- **Ready-to-use** – Collection Swabs and Controls Included
- **Exceptional Customer Service** – Online training modules and technical support staffed by certified medical technologists/professionals

#### Performance Characteristics


	Sensitivity*	Specificity*	Sample Type	Limit of Detection**	Storage Temperature	Time to Result
FLU A	90.3%	96.7%	Nasal/NP Swab	7.1x10 <sup>1</sup> TCID <sub>50</sub> /mL	15-30°C/59-86°F	10 minutes
FLU B	88.0%	99.2%				



#### SENSITIVITY




#### SPECIFICITY




# OSOM® Ultra Plus Flu A&B Test

**SIMPLE TO USE / RESULTS IN MINUTES**



Kit includes one Influenza A+ Control Swab and one Influenza B+ Control Swab

## INTERPRETATION OF RESULTS



## Ordering Information

	Configuration	Catalog Number	CPT Codes
<b>OSOM ULTRA PLUS FLU A&amp;B TEST KIT</b>	25 tests per kit	1032	Flu A: 87804QW
<b>OSOM ULTRA PLUS FLU A&amp;B CONTROL KIT</b>	5 Flu A Swabs 5 Flu B Swabs	1034	Flu B: 87804QW-59

To learn more about all of our OSOM Tests scan the code or go to [osomtests.com](http://osomtests.com)



### THE AMERICAS

SEKISUI Diagnostics, LLC  
One Wall Street  
Burlington, MA 01803  
Phone: 800 332 1042  
Fax: 800 762 6311  
Email: [questions@sekisui-dx.com](mailto:questions@sekisui-dx.com)

\*When compared with RT-PCR

\*\*LOD for Influenza A/Michigan/45/15 (H1N1). Refer to IFU for additional LODs.

© 2021 SEKISUI Diagnostics, LLC. All rights reserved. OSOM® is a registered trademark of SEKISUI Diagnostics, LLC. Because every result matters™ is a trademark of Sekisui Diagnostics, LLC.



**SEKISUI**  
DIAGNOSTICS

Because every result matters™

[sekisuidiagnostics.com](http://sekisuidiagnostics.com)

80-8735-00-00 09/21