

## CRITERION™ EARLOOP MASKS

# RELIABILITY MEETS COMFORT

ASTM RATED PROTECTION IN A MASK  
DESIGNED FOR MAXIMUM WEARABILITY



STYLISH BLACK MASKS NOW AVAILABLE  
IN ASTM LEVEL 1 AND LEVEL 3



For more information, see product Instructions for Use.  
Order today at [HENRYSCHEINDENTAL.COM/CRITERION](https://HENRYSCHEINDENTAL.COM/CRITERION)

**NOW IN BLACK!\***

## CRITERION™ EARLOOP MASKS LEVELS 1, 2 & 3

- **Quality:** ASTM-rated, single use, three-ply masks.
- **Secure Fit:** Ear loops and flexible nosepiece.
- **Gentle On Skin:** Soft, nonirritating inner and outer layer. Not made with natural rubber latex.
- **Comfortable:** Pleated design for enhanced breathability.



\*L1 and L3 only

### ASTM Chart

<b>LOW BARRIER PROTECTION</b>	<b>MODERATE BARRIER PROTECTION</b>	<b>MAXIMUM BARRIER PROTECTION</b>
For general use. Not used for aerosols, spray, or fluids.	Low-mid levels of aerosols, spray, fluids.	High levels of aerosols, spray, fluids.

ASTM REQUIREMENTS	LEVEL 1	LEVEL 2	LEVEL 3
Particle Filtration Efficiency (PFE)	≥ 95% @ 0.1 micron	≥ 98% @ 0.1 micron	≥ 98% @ 0.1 micron
Bacteria Filtration Efficiency (BFE)	≥ 95% @ 3.0 microns	≥ 98% @ 3.0 microns	≥ 98% @ 3.0 microns
Breathability (Delta P)	<5 mm H <sub>2</sub> O/cm <sup>2</sup>	<6 mm H <sub>2</sub> O/cm <sup>2</sup>	<6 mm H <sub>2</sub> O/cm <sup>2</sup>
Fluid Resistance	≥ 80 mm Hg	≥ 120 mm Hg	≥ 160 mm Hg
Flame Spread	Class 1	Class 1	Class 1

## CRITERION™ EARLOOP MASKS 50/BOX

LEVEL 1		LEVEL 2		LEVEL 3	
Blue	104-3809	Blue (US Made)	104-7799	Blue	570-2044
Blue (US Made)	570-2443	Lavender (US Made)	104-8071	Blue (US Made)	570-2448
Lavender	104-8600	Pink (US Made)	104-9200	Lavender	570-2047
Pink	104-3730			Pink	570-2045
Yellow	104-2849			Yellow	570-2046
White	104-6611			White	570-2043
Black <b>NEW!</b>	570-3272			Black <b>NEW!</b>	570-3273
				Teal	570-2042



**BLUE**      **BLUE (US MADE)**      **LAVENDER**      **PINK**      **YELLOW**      **WHITE**      **BLACK**      **TEAL**

Same great product, new packaging! Customers may receive current and new packaging during this transition.

Order today at [HENRYSCHEINDENTAL.COM/CRITERION](https://HENRYSCHEINDENTAL.COM/CRITERION)