

## CRITERION™ NITRILE EXAM N300 GLOVES

# ROBUST PROTECTION IN STYLISH BLACK



**PLUS, NOW FENTANYL-TESTED**  
Exclusive to New Nitrile Exam N200 & N300 Black Gloves

(SEE BACK FOR MORE INFORMATION)



For more information about the full line of Criterion Personal Protective Equipment, scan here or visit [HENRYSCHEINDENTAL.COM/CRITERION](https://HENRYSCHEINDENTAL.COM/CRITERION)

## CRITERION™ NITRILE EXAM N300 GLOVES

- Ultrathin formulation enhances tactile sensitivity and dexterity
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Beaded cuff. 2.0-mil/0.05mm thickness at palm
- Eco-friendly packaging, with 300 gloves in a box, reduces waste and saves space\*



**NEW!**  
**NOW IN BLACK**

\*Excludes X-Large Gloves (570-3486)

### CHEMOTHERAPY TESTED TO ASTM D6978 STANDARD

**Caution:** gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of chemicals used. This glove has been tested for resistance to permeation of various chemotherapy drugs per ASTM D6978, "standard practice for assessment of resistance of medical gloves to permeation by chemotherapy drugs."

**Warning: Do not use Carmustine (BCNU) (3.3mg/ml) and Thiotepa (10mg/ml)** When chemotherapy drugs are present, gloves selection should be based on the specific types of drugs used. Users should review drug labeling or material safety data sheets for the drugs being used to determine an adequate level of protection.

### FENTANYL-TESTED

Tested using ASTM D6978 test method, which evaluates glove permeation resistance to fentanyl. The results show that gloves have resistance to fentanyl for at least 240 minutes.

**Fentanyl breakthrough time (100 mcg/2ml): >240 minutes**

### CHEMOTHERAPY DRUGS & CONCENTRATION

CHEMOTHERAPY DRUGS & CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME
Cisplatin (1.0mg/ml)	>240min
Cyclophosphamide (Cytoxan) (20mg/ml)	>240min
Dacarbazine (DTIC) (10.0mg/ml)	>240min
Doxorubicin Hydrochloride (2.0mg/ml)	>240min
Etoposide (Toposar) (20.0mg/ml)	>240min
Fluorouracil (50.0mg/ml)	>240min
Methotrexate (25mg/ml)	>240min
Mitomycin C (0.5mg/ml)	>240min
Paclitaxel (Taxol) (6.0mg/ml)	>240min
Vincristine Sulfate (1.0mg/ml)	>240min

## BLACK GLOVES

X-Small	300/box	570-3482	Large	300/box	570-3485
Small	300/box	570-3483	X-Large	250/box	570-3486
Medium	300/box	570-3484			

## SHOP THE FULL HEAD-TO-TOE PROTECTION IN STYLISH BLACK



Temple Grip Protective Eyewear



Earloop Masks Levels 1 & 3



Nitrile Exam N200 & N300 Gloves



SMS Jackets and Lab Coats

